

STERIS

INSTRUMENT REPROCESSING | OPERATING THEATRE | ENDOSCOPE CARE



TABLE OF CONTENTS

	000 - Procedural Support
	014 - Instrument Cleaning and Maintenance
	064 - Wraps and Liners
	080 - Instrument Protection
	096 - Tags
	128 - Trays and Accessories
	154 - Container Accessories
	182 - Instrument Identification
	228 - Quality Monitoring and Assurance
	234 - Documentation
	262 - CSSD Equipment

PLACING ORDERS AND PRODUCT ADVICE

If you require advice or wish to place an order, please contact your STERIS representative.

PROCEDURAL SUPPORT

Our products for Procedural Support are designed to help improve a wide range of procedures in the theatre. Our newest product for electrosurgery procedures is Electro Lube. During electrocauterization, eschar build up is common and constant removal of the buildup and handbacks can result in longer procedure times as well as reduced performance of the instrumentation. Electro Lube is anti-stick solution that is used to coat the working end of the instrument which helps keep build up at bay and the instrument from charring.

Anti-Fog is a perfect solution to help reduce condensation and fog on your laparoscopic and arthroscopic lens. Other products such as skin markers, cautery tip cleaner, burr brushes, magnetic mat and drapes, and scrub brushes are also available.

Electro Lube

Anti-Fog

Cautery Tip Cleaner

Burr Brushes

Magnetic Mat and Drapes

Scrub Brushes



1 dispense pouch contents onto sterile field




2 remove paper backing from foam pad and place foam pad onto sterile drape. Note: The foam pad includes a barium filament for radiopacity



3 remove cap and apply Electro Lube onto foam pad or dip electrode tip into bottle



4 wipe electrode tip on foam pad to leave a thin film. Repeat as needed.



ELECTRO LUBE®

An anti-stick solution designed to be used on electrosurgical instrumentation to prevent tissue sticking during the procedure. Also helps protect the instrument from charring. Sterile. Single use.

☆ A

REF	DESCRIPTION	QTY
EL101	4ml bottle, foam pad	20 pcs

A



☆ Top selling items

ANTI-FOG®

Don't let fog affect your sight during a laparoscopic procedure. Anti-Fog is designed to help reduce condensation and fog on your laparoscopic and arthroscopic lens. The answer to a fog-free procedure has never been clearer.

REF	DESCRIPTION	QTY
AF200	6ml anti-fog solution, incl. 50 x 50mm foam pad	20 pcs

☆ A



☆ Top selling items



MAGNETIC INSTRUMENT MAT AND DRAPES

Our reusable magnetic mat provides a safe transfer zone for surgical sharps during procedures. Fluid-resistant material with raised magnetic pods holds stainless steel instruments in place during transfer between staff. Packaged non-sterile.

Our single-use magnetic drapes provide a safe transfer zone for surgical sharps during procedures. Fluid-resistant material with raised magnetic pods holds stainless steel instruments in place during transfer between staff. Packaged sterile.



A

REF	DESCRIPTION	COLOR	SIZE	QTY
MG300400R	reusable magnetic mat	green	300 x 400mm	3 pcs
MG250400	single use magnetic drape	blue	250 x 400mm	30 pcs
MG115076	single use magnetic drape	blue	115 x 76mm	50 pcs

B

C



C



☆ Top selling items



BURR CLEANING BRUSHES

Specially designed for cleaning burs, rasps, files, and saws. Sturdy plastic handle with stainless steel wire bristles. Non-sterile. Single use.

	REF	DESCRIPTION	SIZE	QTY
A	BR101050	round burr, stainless steel bristles	63.5mm	50 pc
B	BR103150	single row, stainless steel bristles	165mm	50 pc
☆ C	BR103350	three rows, stainless steel bristles	165mm	50 pc

CAUTERY TIP CLEANER

Assists in removing material from electrosurgical blades, pencil tips and/or mono-polar, and bi-polar cautery probes during surgery. Foam pad with abrasive surface features adhesive backing for secure, universal placement. Radiopaque. Packaged sterile. Single use.

	REF	DESCRIPTION	COLOR	WIDTH	LENGTH	QTY
☆ D	TC200	foam pad with black abrasive surface	blue	50mm	50mm	100 pcs

☆ Top selling items



HAND AND NAIL BRUSHES

Scrub brushes with a molded block handle to provide secure grip when wet.

	REF	DESCRIPTION	QTY
A	BR4300	medium textured bristles	12 pcs
B	BR4310	hard textured bristles	12 pcs
C	BR4320	medium textured bristles	12 pcs
	BR4200	fine and stiff bristles	1 pcs
☆ D	09391	scrub brush with sponge	40 pcs

SAFETY PINS

Our safety pins are made of stainless steel and are available in lengths of 25, 38, and 50mm. They are packaged into pairs. Sterile.

	REF	LENGTH	QTY
E	06007	25mm	100 pcs
	06008	38mm	100 pcs
	06009	50mm	100 pcs



☆ Top selling items

INSTRUMENT CLEANING AND MAINTENANCE

Our Instrument Cleaning line is designed to provide you with the tools needed to clean your surgical instruments.

With bristles available in nylon or stainless steel, we also provide an incredible variety of lengths and diameters to fit your needs.

For organization of your cleaning brushes, our brush organizer is perfect for hanging in the department. Our stainless steel organizer is easy to keep clean and helps your team choose the right brush in a more efficient manner.

With our vast assortment of brushes - from channel cleaning to endoscope cleaning to sterilizer cleaning - we can help you match the right brush to the right cleaning task at hand.

Toothbrush-Style Brushes

Sterilizer Cleaning Brushes

Channel Cleaning Brushes

Pipe Cleaners

Endoscope Brush

Preparation Aids





TOOTHBRUSH-STYLE CLEANING BRUSHES

Our toothbrush-style brushes are ideal for cleaning the outside of the surfaces of instruments before automated disinfection and sterilization. Regardless of the stubbornness of the soiling of your instruments to be cleaned, you are sure to find the right solution thanks to the varying degrees of hardness, sizes, and shapes of the brushes. Even coagulated blood residues can be easily removed due to specially developed nylon bristles. Our specialized brushes are also suitable for cleaning pan burrs and removing residues in the joints of needle holders and scissors.

TOOTHBRUSH-STYLE CLEANING BRUSHES

	REF	DESCRIPTION	BRUSH HEAD LENGTH	BRISTLE LENGTH	WORKING LENGTH	QTY
A	09320	plastic bristles for stubborn contamination - extra hard	30/45mm	12/12mm	165mm	10 pcs
B	09325	stainless steel bristles for stubborn contamination - extra hard	40mm	15mm	180mm	1 pc
☆ C	S210050	stainless steel, three rows				50 pcs
	09459	plastic bristles for stubborn contamination - extra hard	40/30mm	10/10mm	175mm	2 pcs
	N222150	extra rigid nylon, double-ended, three rows/one row				50 pcs
D	09478	plastic bristles - soft	75mm	15mm	215mm	1 pc
	09485	plastic bristles for stubborn contamination, angled handle - hard	75mm	15mm	235mm	2 pcs
☆ E	N2200	extra rigid nylon, angled handle				2 pcs
☆ F	26491	plastic bristles - medium soft	75mm	15mm	215mm	1 pc
G	26492	plastic bristles - hard	75mm	15mm	215mm	1 pc
☆ H	09488	plastic bristles	40mm	15mm	180mm	1 pc
	N210050	nylon, three rows	-	-	-	50 pcs

☆ Top selling items



A



B



C



D



E



F

TOOTHBRUSH-STYLE CLEANING BRUSHES

Brushes with nylon or stainless steel bristles on a sturdy, plastic handle. Designed for external cleaning.

	REF	DESCRIPTION	BRISTLE LENGTH	WORKING LENGTH	QTY
A	N3000	made of nylon	9.5mm	178mm	1 pc
B	N4000	made of nylon	9.5mm	184mm	1 pc
C	N212150	nylon, double-ended, three rows/one row	10/5mm	175mm	50 pcs
D	09098	hard nylon bristles	10/5mm	175mm	2 pcs
E	S212150	stainless steel, double-ended, three rows/one row	40	177mm	50 pcs

CLEANING BRUSH

	REF	DESCRIPTION	BRISTLE LENGTH	BRUSH HEAD LENGTH	WORKING LENGTH	QTY
E	BR5450	brush bristles made of nylon, stainless steel wire	13/25mm	50/82mm	300mm	1 pc

☆ Top selling items



A



B



C




D

STERILIZER CLEANING BRUSHES

Nylon bristles on twisted stainless steel handle with rosebud tip or on a sturdy plastic handle.


	REF	DESCRIPTION	FORMAT Ø	LENGTH	QTY
A	BR5455	finger loop brush, made of nylon	25 x 270mm	267mm	1 pc
B	BR5500	sterilizer cleaning brush, made of nylon with a plastic handle	35 x 710mm	710mm	1 pc
C	BR4000	with nylon bristles, plastic handle	41 x 200mm	410mm	1 pc
D	BR5000	sterilizer cleaning brush, made of polyester	152 x 127mm	800mm	1 pc



CHANNEL CLEANING BRUSHES

Nylon bristles, twisted stainless steel wire.

	REF	COLOR	Ø DIA	BRUSH HEAD LENGTH	WORKING LENGTH	QTY
A	09095	-	1mm	13mm	450mm	2 pcs
B	09319	-	1.8mm	13mm	450mm	2 pcs
C	45906	white	1.6mm	22mm	215mm	3 pcs
D	45908	black	2mm	21mm	215mm	3 pcs
E	45910	bright pink	2.5mm	21mm	215mm	3 pcs
F	45912	orange	3mm	18mm	215mm	3 pcs
G	45914	blue	4mm	18,5mm	215mm	3 pcs
H	BR0808750		2.2mm	13mm	200mm	50 pcs




BRUSH TIP DETAILS SHOWN ACTUAL SIZE

CHANNEL CLEANING BRUSHES

Channel cleaning brushes featuring nylon bristles on a twisted stainless steel handle with loop end.

	REF	Ø DIA	BRUSH HEAD LENGTH	WORKING LENGTH	QTY
A	BR0603950	1mm	6.35mm	152mm	50 pcs
B	BR0803950	1mm	12.7mm	203.2mm	50 pcs
	BR1803950	1mm	6.35mm	457mm	50 pcs
C	BR0606250	1.58mm	12.7mm	152mm	50 pcs
	BR1806250	1.58mm	12.7mm	457mm	50 pcs
D	BR0807050	1.77mm	12.7mm	203.2mm	50 pcs
	BR1807050	1.77mm	12.7mm	457mm	50 pcs
	BR2407050	1.77mm	12.7mm	609.6mm	50 pcs
E	BR0607950	2mm	12.7mm	152mm	50 pcs
F	BR0807950	2mm	19.05mm	203.2mm	50 pcs
	BR1207950	2mm	19.05mm	305mm	50 pcs
G	BR2408050	2mm	25.4mm	609.6mm	50 pcs
H	BR0808750	2.2mm	12.7mm	203.2mm	50 pcs
	BR1808750	2.2mm	12.7mm	457mm	50 pcs
I	BR0809250	2.3mm	12.7mm	203.2mm	50 pcs
J	BR2409850	2.49mm	12.7mm	609.6mm	50 pcs
K	BR0810550	2.67mm	25.4mm	203.2mm	50 pcs
	BR1810550	2.67mm	25.4mm	457mm	50 pcs
L	BR0811850	3mm	38.1mm	203.2mm	50 pcs
M	BR1211850	3mm	69.85mm	305mm	50 pcs



BRUSH TIP DETAILS SHOWN ACTUAL SIZE

A

B

C

D

E

F

G

CHANNEL CLEANING BRUSHES

Channel cleaning brushes featuring nylon bristles on a twisted stainless steel handle with loop end.



A

B

C

D




E

F

G

REF	Ø DIA	BRUSH HEAD LENGTH	WORKING LENGTH	QTY
BR1812550	3.18mm	50.8mm	457mm	50 pcs
BR0813150	3.3mm	19.05mm	203.2mm	50 pcs
BR0814450	3.66mm	19.05mm	203.2mm	50 pcs
BR0815850	4mm	19.05mm	203.2mm	50 pcs
BR1215850	4mm	19.05mm	305mm	50 pcs
BR1815850	4mm	19.05mm	457mm	50 pcs
BR1219750	5mm	76.2mm	305mm	50 pcs
BR1619750	5mm	76.2mm	406mm	50 pcs
BR1819750	5mm	76.2mm	457mm	50 pcs
BR2419750	5mm	76.2mm	609.6mm	50 pcs
BR2420050	5.08mm	50.8mm	609.6mm	50 pcs
BR1223650	6mm	50.8mm	305mm	50 pcs
BR1623650	6mm	50.8mm	406mm	50 pcs
BR2423650	6mm	50.8mm	609.6mm	50 pcs

☆ Top selling items



BRUSH TIP DETAILS SHOWN ACTUAL SIZE

A

B

C

D

E

F

G

CHANNEL CLEANING BRUSHES

Channel cleaning brushes featuring nylon bristles on a twisted stainless steel handle with loop end.

	REF	Ø DIA	BRUSH HEAD LENGTH	WORKING LENGTH	QTY
★ A	BR0625050	6.35mm	50.8mm	152mm	50 pcs
	BR0825050	6.35mm	50.8mm	203.2mm	50 pcs
★ B	BR1825050	6.35mm	50.8mm	457mm	50 pcs
	BR1627650	7mm	76.2mm	406mm	50 pcs
★ C	BR2427650	7mm	76.2mm	609.6mm	50 pcs
	BR2431250	7.9mm	76.2mm	609.6mm	50 pcs
★ D	BR0637550	9.5mm	38.1mm	152mm	50 pcs
★ E	BR2437550	9.5mm	76.2mm	609.6mm	50 pcs
★ F	BR1639350	10mm	76.2mm	406mm	50 pcs
	BR1839350	10mm	76.2mm	457mm	50 pcs
	BR2439350	10mm	76.2mm	609.6mm	50 pcs
★ G	BR2442550	10.8mm	76.2mm	609.6mm	50 pcs

★ Top selling items



A

B

C

D


E

BRUSH TIP DETAILS SHOWN ACTUAL SIZE

CHANNEL CLEANING BRUSHES

Channel cleaning brushes featuring nylon bristles on a twisted stainless steel handle with loop end.

REF	Ø DIA	BRUSH HEAD LENGTH	WORKING LENGTH	QTY
BR0844050	11.18mm	76.2mm	203.2mm	50 pcs
BR2444050	11.18mm	76.2mm	609.6mm	50 pcs
BR1644350	11.25mm	76.2mm	406mm	50 pcs
BR2444350	11.25mm	76.2mm	609.6mm	50 pcs
BR1647250	12mm	76.2mm	406mm	50 pcs
BR0650050	12.7mm	50.8mm	152mm	50 pcs
BR1659150	15mm	50.8mm	406mm	50 pcs
BR1859150	15mm	50.8mm	457mm	50 pcs



BRUSH TIP DETAILS SHOWN ACTUAL SIZE

Instrument Cleaning and Maintenance

Channel Cleaning Brushes

CHANNEL CLEANING BRUSHES


Channel cleaning brushes featuring nylon bristles on a twisted stainless steel handle with loop end.

	REF	Ø DIA	BRUSH HEAD LENGTH	WORKING LENGTH	QTY
A	BR0863050	16mm	50.8mm	203mm	50 pcs
B	BR1678750	20mm	50.8mm	406mm	50 pcs
C	BR1888050	22.35mm	88.9mm	457mm	50 pcs
D	BR2488050	22.35mm	107.95mm	609.6mm	50 pcs


CHANNEL CLEANING BRUSHES

Channel cleaning brushes featuring nylon bristles on a twisted stainless steel handle with loop end (pull through style excluded).

REF	DESCRIPTION	Ø DIA	BRUSH HEAD LENGTH	WORKING LENGTH	QTY
BR3044	nylon	0.5mm	6.35mm	76mm	10 pcs
BR3045	nylon	0.5mm	6.35mm	203.2mm	10 pcs
BR3047	nylon	0.9mm	12.7mm	101.6mm	10 pcs
BR412050	nylon, centered (pull through)	6.35 to 15.24mm	50.8mm	457.2mm	50 pcs
BR4140	nylon, centered (pull through)	12.7mm	76.2mm	609.6mm	6 pcs



BRUSH TIP DETAILS SHOWN ACTUAL SIZE




BRUSH TIP DETAILS SHOWN ACTUAL SIZE

CHANNEL CLEANING BRUSHES

Channel cleaning brushes featuring nylon bristles on a twisted stainless steel handle with loop end.

	REF	Ø DIA	BRUSH HEAD LENGTH	WORKING LENGTH	QTY
A	BR0811850	3mm	38mm	200mm	50 pcs
B	BR1812550	3mm	51mm	450mm	50 pcs
C	BR0814450	3.7mm	19mm	203mm	50 pcs
D	BR1815850	4mm	19mm	450mm	50 pcs
E	BR2419750	5mm	76mm	610mm	50 pcs
F	BR1888050	22mm	108mm	450mm	50 pcs
G	BR2488050	22mm	108mm	610mm	50 pcs



BRUSH TIP DETAILS SHOWN ACTUAL SIZE

CHANNEL CLEANING BRUSHES

Channel cleaning brushes featuring nylon bristles on a twisted stainless steel handle with loop end.

	REF	Ø DIA	BRUSH HEAD LENGTH	WORKING LENGTH	QTY
A	BR811050	9.5mm	25mm	95mm	50 pcs
B	BR8115	13mm	101mm	203mm	3 pcs
C	BR821050	4.8mm	12,7mm	178mm	50 pcs
D	BR8410	13mm	76,2mm	152mm	3 pcs



A



B



C



D

BRUSH TIP DETAILS SHOWN ACTUAL SIZE

CHANNEL CLEANING BRUSHES


Nylon bristles, twisted stainless steel wire and protective tip.

	REF	DESCRIPTION	Ø DIA	BRUSH HEAD LENGTH	WORKING LENGTH	QTY
A	09303	nozzle brush, made of nylon	10mm	45mm	140mm	5 pcs
B	09490	cleaning brush for medullary reamer brush tip deburred, V2A wire	3mm	100mm	800mm	3 pcs
C	09491	for medullary reamers, V2A stainless wire tip with cap	7mm	100mm	800mm	3 pcs
D	09492	for medullary reamers, stainless steel wire (V2A), tip with end cap	5mm	100mm	800mm	3 pcs


CHANNEL CLEANING BRUSHES

Nylon bristles, twisted stainless steel wire.

REF	DESCRIPTION	Ø DIA	BRUSH HEAD LENGTH	WORKING LENGTH	QTY
09331	for rectoscopes	20mm	150mm	450mm	5 pcs
09366	front 90 mm without trim	0.6mm	8mm	450mm	5 pcs
09388	front 90 mm without trim	1.4mm	8mm	450mm	5 pcs
09338	with plastic handle, not suitable for washer-disinfectors	3.5mm	60mm	265mm	2 pcs



BRUSH TIP DETAILS SHOWN ACTUAL SIZE



BRUSH TIP DETAILS SHOWN ACTUAL SIZE

CLEANING BRUSHES FOR MINIMALLY INVASIVE SURGICAL INSTRUMENTS

Nylon bristles, twisted stainless steel wire with protective tip.

REF	DESCRIPTION	Ø DIA	BRUSH LENGTH	WORKING LENGTH	QTY
A 09050	nozzle brush, made of nylon	2.5mm	30mm	500mm	5 pcs
B 09051	nozzle brush, made of nylon	7mm	80mm	350mm	5 pcs
☆ C 09052	nozzle brush, made of nylon	4mm	100mm	300mm	5 pcs
☆ D 09053	nozzle brush, made of nylon	6mm	100mm	300mm	5 pcs
☆ E 09054	nozzle brush, made of nylon	8mm	100mm	300mm	5 pcs
F 09055	nozzle brush, made of nylon	10mm	100mm	300mm	5 pcs
☆ G 09056	nozzle brush, made of nylon	15mm	100mm	300mm	5 pcs
☆ H 09059	nozzle brush, made of nylon	30mm	80mm	300mm	5 pcs
I 09058	nozzle brush, made of nylon	7mm	50mm	405mm	5 pcs


GENERAL INSTRUMENT CLEANING BRUSH

With nylon bristles, non heat-resistant, not suitable for washer-disinfectors.

REF	DESCRIPTION	Ø DIA	BRUSH LENGTH	WORKING LENGTH	QTY
J 09064	double-sided	10/15mm	30/12-15mm	155mm	10 pcs



☆ Top selling items



BRUSH TIP DETAILS SHOWN ACTUAL SIZE

- A
- B
- C
- D
- E
- F
- G
- H
- I
- J
- K
- L


CLEANING BRUSHES FOR MINIMALLY INVASIVE SURGICAL INSTRUMENTS

Nylon bristles, twisted stainless steel wire with protective tip (where applicable).

- A
- ☆ B
- C
- D
- ☆ E
- F
- G
- ☆ H
- I
- ☆ J
- K
- ☆ L

REF	Ø DIA	BRUSH LENGTH	WORKING LENGTH	QTY
09308	2.5mm	50mm	300	5 pcs
09313	2.5mm	50mm	610	5 pcs
09306	5mm	100mm	300	5 pcs
09312	5mm	100mm	610	5 pcs
09311	7mm	100mm	610	5 pcs
09316	7mm	100mm	650	5 pcs
09068	10mm	100mm	500	5 pcs
09310	10mm	100mm	610	5 pcs
09315	10mm	100mm	650	5 pcs
09309	12mm	100mm	610	5 pcs
09314	12mm	100mm	650	5 pcs
09083	50mm	100mm	1400	5 pcs

☆ Top selling items



BRUSH TIP DETAILS SHOWN ACTUAL SIZE

CLEANING BRUSHES FOR MINIMALLY INVASIVE SURGICAL INSTRUMENTS

Nylon bristles, twisted stainless steel wire with protective tip (where applicable).

	REF	Ø DIA	BRUSH LENGTH	WORKING LENGTH	QTY
A	09339	2mm	100mm	300mm	5 pcs
B	09342	2mm	60mm	350mm	5 pcs
C	09340	2.5mm	100mm	300mm	5 pcs
D	09343	2.5mm	80mm	400mm	5 pcs
E	09341	3mm	100mm	300mm	5 pcs
F	09344	3mm	100mm	450mm	5 pcs
G	09335	3.5mm	100mm	250mm	5 pcs
H	09346	3.5mm	35mm	500mm	5 pcs
I	09345	4mm	100mm	500mm	5 pcs
J	09317	5mm	100mm	650mm	5 pcs
K	09318	5mm	80mm	405mm	5 pcs
L	09324	5mm	35mm	500mm	5 pcs
M	09336	5mm	100mm	250mm	5 pcs
N	09332	20mm	100mm	500mm	5 pcs
O	09334	40mm	60mm	500mm	5 pcs



CLEANING BRUSHES WITH FAN TIP

Channel cleaning brushes featuring a fan tip to assist with cleaning closed-end channels. Nylon bristles on a twisted stainless steel handle with loop end. Also available as a pack of 50.

	REF	Ø DIA	BRUSH HEAD LENGTH	WORKING LENGTH	QTY
A	FT12197	5mm	50mm	300mm	2 pcs
A	FT18197	5mm	50mm	460mm	2 pcs
A	FT24197	5mm	50mm	610mm	2 pcs
B	FT12250	6.4mm	50mm	300mm	2 pcs
B	FT18250	6.4mm	50mm	460mm	2 pcs
B	FT24250	6.4mm	50mm	460mm	2 pcs
C	FT12276	7mm	50mm	460mm	2 pcs
C	FT18276	7mm	50mm	460mm	2 pcs
C	FT24276	7mm	50mm	460mm	2 pcs
D	FT12375	9.5mm	76mm	460mm	2 pcs
D	FT18375	9.5mm	76mm	460mm	2 pcs
D	FT24375	9.5mm	76mm	460mm	2 pcs
E	FT12440	11mm	76mm	460mm	2 pcs
E	FT18440	11mm	76mm	460mm	2 pcs
E	FT24440	11mm	76mm	460mm	2 pcs

BONE REAMER BRUSHES

Channel cleaning brushes featuring nylon bristles on a rigid stainless steel handle with loop end.


	REF	Ø DIA	BRUSH HEAD LENGTH	WORKING LENGTH	QTY
F	BR31118	3mm	76mm	780mm	3 pcs
G	BR31197	5mm	76mm	780mm	3 pcs
H	BR31276	7mm	76mm	780mm	3 pcs

ACETABULAR REAMER BRUSH

Designed for cleaning acetabular bone reamer heads. Rigid with nylon bristles. Double-ended. Round pigtail style brush on one end, straight precision style brush on the other.

	REF	Ø DIA	BRUSH HEAD LENGTH	WORKING LENGTH	QTY
☆ I	N20003	40mm	45mm	190mm	3 pcs

☆ Top selling items




CLEANING BRUSHES FOR INSTRUMENT CAVITIES

Nylon bristles on a twisted stainless steel handle with loop end (pull through style excluded), and single use plastic shaft brushes.

	REF	Ø DIA	SHAFT DIAMETER	WORKING LENGTH	QTY
A	26475	5 x 50mm	3.2mm	610mm	5 pcs
B	26476	7 x 50mm	3.2mm	610mm	5 pcs
C	26478	15 x 50mm	4mm	610mm	5 pcs
D	BR412050	6 x 15mm	-	45mm	50 pcs
E	BR4140	13mm	-	610mm	6 pcs

PULL THROUGH



ACRYLIC TIP CLEANING BRUSHES

Channel cleaning brushes featuring a protective acrylic tip for an extra level of protection during cleaning. Nylon bristles on a flexible handle. No loop end.

	REF	Ø DIA	BRUSH HEAD LENGTH	WORKING LENGTH	QTY
A	AT12062	1.6mm	25mm	300	2 pcs
	AT18062	1.6mm	25mm	457	2 pcs
	AT24062	1.6mm	25mm	610	2 pcs
B	AT12087	2.2mm	25mm	300	2 pcs
C	AT18087	2.2mm	25mm	457	2 pcs
	AT24087	2.2mm	25mm	610	2 pcs
	AT12105	2.7mm	25mm	300	2 pcs
D	AT18105	2.7mm	25mm	457	2 pcs
	AT24105	2.7mm	25mm	610	2 pcs
E	AT12118	3mm	25mm	300	2 pcs
F	AT18118	3mm	25mm	457	2 pcs
G	AT24118	3mm	25mm	610	2 pcs
	AT12125	3.2mm	25mm	300	2 pcs
H	AT18125	3.2mm	25mm	457	2 pcs
	AT24125	3.2mm	25mm	610	2 pcs
	AT12131	3.3mm	25mm	300	2 pcs
I	AT18131	3.3mm	25mm	457	2 pcs
	AT24131	3.3mm	25mm	610	2 pcs
J	AT12158	4mm	38mm	300	2 pcs
K	AT18158	4mm	38mm	457	2 pcs
	AT24158	4mm	38mm	610	2 pcs
L	AT12197	5mm	51mm	300	2 pcs
M	AT18197	5mm	51mm	457	2 pcs
	AT24197	5mm	55mm	610	2 pcs
N	AT12250	6.4mm	51mm	300	2 pcs
	AT18250	6.4mm	51mm	457	2 pcs
	AT24250	6.4mm	55mm	610	2 pcs
O	AT12276	7mm	51mm	300	2 pcs
P	AT18276	7mm	51mm	457	2 pcs
	AT24276	7mm	55mm	610	2 pcs
Q	AT12375	9.5mm	76mm	300	2 pcs
	AT18375	9.5mm	76mm	457	2 pcs
	AT24375	9.5mm	76mm	610	2 pcs
R	AT12440	11mm	76mm	300	2 pcs
	AT18440	11mm	76mm	457	2 pcs
	AT24440	11mm	76mm	610	2 pcs



- A
- B
- C
- D
- E
- F
- G
- H
- I

- A
- B
- C
- D
- E

Instrument Cleaning and Maintenance

Pipe Cleaners

FLEXISTEM AND MICROCLENE

Flexistem and Microclene are unique cleaners for all small-bore cannula. Various diameters from 3mm to 15mm are available, constructed from very soft to hard bristles. Single use. Flexistem FLEX007A and FLEX007B has a wavy bristle configuration to reduce buckling which can occur in long instruments.


FLEXISTEM

REF	TYPE	MATERIAL	Ø DIA	LENGTH	QTY
FLEX001	soft cleaner	cotton	3mm	1000mm	1 reel
FLEX003	double ended	polyester/nylon	3mm	30mm	100 pcs
FLEX004	fleece	polyester	9mm	30mm	50 pcs
FLEX005	fleece	polyester	12mm	30mm	50 pcs
FLEX006	bristle stem	polyester	15mm	30mm	20 pcs
FLEX007	bristle stem	polyester/nylon	15mm	500mm	1 reel
FLEX007A	bristle stem	polyester/nylon	15mm	500mm	1 reel
FLEX033	fleece	polyester	3mm	30mm	50 pcs

MICROCLENE

REF	TYPE	MATERIAL	Ø DIA	LENGTH	QTY
FLEX021	super fine	polyester/nylon	3mm	750cm	1 reel
FLEX022	medium fine	polyester/nylon	6mm	750cm	1 reel
FLEX023	super brush	polyester/nylon	3mm	750cm	1 reel
FLEX024	medium brush	polyester/nylon	6mm	750cm	1 reel


☆ Top selling items



ENDOSCOPE BRUSH

Use our DuoSwift Combination Squeegee Brush to work smarter, not harder, during flexible endoscope reprocessing. Designed to clean endoscope channel diameters from 2.8mm - 7.0mm, the unique design of DuoSwift features a tapered, nylon-bristled brush on one end and a series of flexible disks (squeegees) on the other. This brush does the heavy lifting for you; the combination of bristle brush and squeegee provides cleaning action that removes more soil in a single pass.


REF	DESCRIPTION	BRISTLE Ø	SQUEEGEE Ø	STUBBY BRUSH	OVERALL LENGTH	QTY
BX00711618	DuoSwift combination squeegee brush	5 – 7 (tapered)mm	2.5 – 5.2mm	n/a	230mm	50/box
BX00711619	DuoSwift combination squeegee brush with valve brush	5 – 7 (tapered)mm	2.5 – 5.2mm	10mm	230mm	50/box





PREPARATION AIDS - FOR MULTIPLE USE

Holders and band clamps to facilitate cleaning of instruments by holding open the instruments.

	REF	DESCRIPTION	COLOR	MATERIAL	DIAMETER	QTY
A	20413	holder for microclip	-	stainless steel	4.8mm	1 pc
B	20409	holder for clamps	-	stainless steel	6mm	5 pcs
C	20429	holder for clamps	-	stainless steel	8.4mm	5 pcs
D	20410	holder for clamps	-	stainless steel	10mm	5 pcs
E	20430	holder for clamps	-	stainless steel	14mm	5 pcs
F	20460	clamp	white	stainless steel with silicone cover	40mm	1 pc
G	20411	clamp	yellow	stainless steel with silicone cover	50mm	1 pc
H	20408	clamp	blue	stainless steel with silicone cover	60mm	1 pc
I	20406	clamp	red	stainless steel with silicone cover	70mm	1 pc
J	20412	clamp	green	stainless steel with silicone cover	80mm	1 pc
K	20407	holder for instruments with a ring handle	-	stainless steel	-	1 pc





TRAY ORGANIZATION MATS

Suitable for washer-disinfectors, silicone or stainless steel designed to provide support for instruments during high level disinfection. Heat resistant 134 °C.

☆ A

REF	DESCRIPTION	OUTER DIMENSIONS	INNER DIMENSIONS	QTY
09020	for 1/1 strainer baskets, frames rigid made of silicone	445 x 230mm	425 x 210mm	1 pc
08847	stainless steel, rigid mesh	445 x 230mm	-	1 pc

PREPARATION ACCESSORIES

☆ C

REF	DESCRIPTION	FORMAT	QTY
20459	Rack, fitting 8 kidney dishes	460 x 220 x 165mm	1 pc

☆ Top selling items

WRAPS AND LINERS

Our wraps and liners are available in a wide variety of sizes but can also be customized to your specific needs.

Barrier Wraps provide a bacterial barrier and protect the contents of a set or pack from contamination. When classed as alcohol repellent, barrier wraps can be used within the sterile field.

Choice Wrap is designed to reduce the risk of rips, tears, and wet packs. Choice SMS are made of polypropylene. Choice Wraps can be interleaved in blue or green to meet the sterilization service user requirements.

All Transport Wraps must be used in conjunction with one compliant Barrier Wrap to completely prevent the ingress of microorganisms. Use of a single barrier layer with Transport Wrap can be a cost-effective solution that reduces damage in handling and transport. If adopted, changes to both wrapping and unwrapping procedures may be required to ensure aseptic opening.

Choice Wrap

Crepe Wrap

Tray Liners

Soaker Sheets

Adhesive Tape and Dispenser

Transport Wrap





A



B



CHOICE WRAP 55gsm

A fully synthetic SMS product designed to offer flexibility and enable users to pack conventionally or with ‘two layers’ at once. Choice Wrap is strong and available with alcohol repellence in both green and blue for use within the sterile field. Wraps are also available with blue and green interleaved. Compatible with steam sterilization. Single use.

A CHOICE WRAP BLUE

REF	FORMAT	QTY
CB060	60 x 60cm	350 pcs
CB075	75 x 75cm	300 pcs
CB090	90 x 90cm	300 pcs
CB100	100 x 100cm	200 pcs
CB120	120 x 120cm	120 pcs
CB1315	130 x 150cm	100 pcs

B CHOICE WRAP GREEN

REF	FORMAT	QTY
CG075	75 x 75cm	300 pcs
CG100	100 x 100cm	200 pcs
CG120	120 x 120cm	120 pcs
CG1315	130 X 150 cm	100 pcs

CHOICE WRAP COMBINATION (GREEN /BLUE)

REF	FORMAT	QTY
CGCB060	60 x 60cm	500 pcs
CGCB075	75 x 75cm	200 pcs
CGCB090	90 x 90cm	200 pcs
CGCB100	100 x 100cm	200 pcs
CGCB120	120 x 120cm	100 pcs
CGCBU075	75 X 75 cm	300 pcs

☆ Top selling items



RELIANCE® CREPE WRAP 60gsm

Crepe wrap is suitable for lighter instruments and sets. Manufactured from 100% medical quality cellulose fibers and combines the most economical wrapping with good drapability. Compatible with steam and EO sterilization. Single use.

☆

A

☆

B

REF	COLOR	FORMAT	QTY
GRNCREPE060N	green	60 x 60cm	750 pcs
GRNCREPE075NU	green	75 x 75cm	174 pcs
GRNCREPE090U	green	90 x 90cm	126 pcs
GRNCREPE100U	green	100 x 100cm	126 pcs
GRNCREPE120U	green	120 x 120cm	126 pcs
WHCREPE050N	white	50 x 50cm	750 pcs
WHCREPE060N	white	60 x 60cm	750 pcs
WHCREPE100U	white	100 x 100cm	252 pcs

A



B



☆ Top selling items

A



A TRAY LINERS

Our Tray Liners have high absorbency and rapid dispersion to remove condensation, encouraging dry packs and improving the efficacy of the sterilization process. To be used above or beneath trays. 80gsm thickness. Compatible with steam sterilization. Single use.



REF	FORMAT	QTY
TLC004	200 x 305mm	2000 pcs
TLC005	250 x 250mm	2400 pcs
TLC006	250 x 500mm	1200 pcs
TLC007	305 x 410mm	1000 pcs
TLC008	457 x 550mm	800 pcs
TLC009	305 x 305mm	2000 pcs
TLC010	305 x 457mm	1000 pcs
TLC011	305 x 610mm	1000 pcs
TLC012	457 x 457mm	1000 pcs
TLC013	457 x 610mm	600 pcs
TLC027	350 x 500mm	1000 pcs
TLC036	305 x 500mm	1000 pcs
TLC045	400 x 500mm	1000 pcs
TLC051	250 x 305mm	2000 pcs
TLC055	250 x 550mm	1200 pcs
TLC085	250 x 480mm	1000 pcs
TLC100	305 x 350mm	1000 pcs
TLC110	150 x 250mm	2000 pcs
TLC340	250 x 340mm	2000 pcs

☆ Top selling items



A SOAKER SHEETS

140gsm Soaker Sheets are for use under trays to protect inner wraps from sharp corners and feet on trays and to absorb condensate during sterilization. Soaker Sheets can also be used as a shelf liner for added protection. Single use.



REF	SIZE	QTY
SOK001SS	750 x 500mm	50 sheets
SOK002SS	500 x 375mm	100 sheets
SOK003SS	610 x 305mm	100 sheets
SOK004SS	500 x 200mm	100 sheets
SOK005SS	500 x 500mm	100 sheets
SOK006SS	250 x 250mm	100 sheets
SOK007SS	500 x 280mm	100 sheets
SOK008SS	1219 x 560mm	50 sheets
SOK009SS	1524 x 560mm	50 sheets

☆ Top selling items



TAPE DISPENSERS

	REF	DESCRIPTION	QTY
A	09001	table dispenser for 1 roll of indicator tape	1 pc
B	09002	table dispenser for 2 rolls of indicator tape	1 pc
	09008	replacement knife for table dispenser REF 09002	1 pc
C	09035	table dispenser for indicator tape up to 25mm strip length	1 pc



A **TRANSPORT COVER BAGS**

Our Transport Cover Bags are made of the same material as the Dry Transport Cover Wraps and can alternatively be sealed with indicator tape. Available in different sizes.

REF	FORMAT	QTY
TCOVER001	400 x 700 x 90mm	200 pcs
TCOVER002	500 x 700 x 90mm	200 pcs
TCOVER003	600 X 700 X 90mm	100 pcs
TCOVER008	550 X 800 X 90mm	100 pcs

B **TRANSPORT WRAP DRY 60gsm**

Protection against damaged packs, with strength and tear resistance with a patented absorbent layer. Particularly suitable for use when handling heavy orthopaedic sets. Available in purple, to distinguish from barrier wraps and also helps to support identification during transport and storage. Compatible with steam sterilization. Single use.

REF	FORMAT	QTY
TDRY090	90 x 90	300 pcs
TDRY120	120 x 120	120 pcs
TWRAP100	100 x100	200 pcs
TWRAP120	120 x 120	120 pcs

B



☆ Top selling items

INSTRUMENT PROTECTION

The protection of instruments has never been more important. As a result of increasing centralization, the merger of CSSDs and external reprocessing, transport routes are becoming ever longer. This means that sensitive and complex instruments have to be protected not only within your own hospital, but often when being transported off-site.

Instrument Sleeves are made of special cardboard or transparent film and are used to protect high-quality instruments from damage during preparation, transport and storage. The shape of the protectors is not directly adapted to the respective instrument. This has the advantage that these protective sleeves are suitable for the most common instruments as well as for special shapes (such as hooks or knives).

Instrument Tip Caps are designed to extend the life of your instruments. Our tip caps protect your valuable instruments from unnecessary wear or damage, helping reduce repair costs while giving you the peace of mind that your instrument investment is being protected. In addition to protecting and extending the life of your valuable instruments, tip caps also help to protect sterilization packs from puncture and protect employees from sharp instruments.

Available in a variety of shapes and sizes to fit all your instrument protection needs, we have the perfect tip cap for each of your various instruments, from flat caps to twin caps to instrument holders, and they are available with or without perforations, depending on your preference. Our tip caps are radiopaque, single use, and compatible with steam sterilization.

Instrument Sleeves

Instrument Protection Caps Round

Instrument Protection Tapered Caps

Instrument Protection Caps Flat

Instrument Retractor Caps

Instrument Protection Twin Caps

Instrument Organizer







INSTRUMENT SLEEVES

Protection sleeves are used to protect high-quality instruments from damage during preparation, transport, and storage. They are therefore suitable for the most common instruments such as hooks or knives. Single use.

- ☆ A
- ☆ B
- ☆ C
- ☆ D

REF	DESCRIPTION	MATERIAL	FORMAT*	QTY
01455	hybrid instrument sleeve	cardboard/film	120 x 32mm	100 pcs
01457	hybrid instrument sleeve	cardboard/film	160 x 40mm	100 pcs
01476	instrument sleeve	film	60 x 10mm	1000 pcs
01471	folding borders can be folded inwards to prevent narrow instruments from slipping out	cardboard	105 x 25mm	1000 pcs
01477	folding borders can be folded inwards to prevent narrow instruments from slipping out	film	140 x 50mm	100 pcs
01479	for amputation knives, with 8 holes	cardboard	123 x 33mm	100 pcs
01478	for amputation knives, with 8 holes	cardboard	246 x 33mm	100 pcs
PPT001	white card		160 x 50mm	2000 pcs
PPT002	white card		160 x 50mm	1000 pcs
PPT004	white card		160 x 50mm	1000 pcs

* Useful length x max. instrument width



☆ Top selling items



PRODUCTS SHOWN ACTUAL SIZE

Guards and caps are used to protect premium instruments from damage during transportation, handling and sterilization. Furthermore, while they don't provide a sterile barrier, they can help to reduce costs due to rips or punctures in sterilization packaging. Standard vinyl protective tip guards are also available in a range of shapes.

ROUND CAPS

For single use.

	REF	DESCRIPTION	COLOR	Ø DIA	QTY
☆ A	30401	perforated	clear	1.6 x 19mm	100 pcs
☆ B	30801	perforated	white transparent	1.6 x 19mm	100 pcs
C	31501	perforated	white	1.6 x 19mm	100 pcs
D	30402	perforated	clear	2.0 x 19mm	100 pcs
☆ E	30802	perforated	blue transparent	2.0 x 19mm	100 pcs
F	31502	perforated	blue	2.0 x 19mm	100 pcs
G	30403	perforated	clear	2.8 x 19mm	100 pcs
☆ H	30803	perforated	green transparent	2.8 x 19mm	100 pcs
I	31503	perforated	green	2.8 x 19mm	100 pcs
J	30404	perforated	clear	3.2 x 25mm	100 pcs
☆ K	30804	perforated	red transparent	3.2 x 25mm	100 pcs
L	31504	perforated	red	3.2 x 25mm	100 pcs
M	30499	perforated	clear	4.7 x 38mm	100 pcs
N	30405	perforated	clear	5.0 x 25mm	100 pcs
☆ O	30805	perforated	yellow transparent	5.0 x 25mm	100 pcs
P	31505	perforated	yellow	5.0 x 25mm	100 pcs
Q	30406	perforated	clear	10 x 19mm	100 pcs
R	30806	perforated	purple transparent	10 x 19mm	100 pcs
S	30703	for MIS instruments	green	4.3 x 29mm	50 pcs
T	30705	for MIS instruments	yellow	4.3 x 51mm	50 pcs
U	01858	for obturators	clear	5.0 x 43mm	20 pcs
V	30709	for MIS instruments	clear	66 x 9.0mm	50 pcs

☆ Top selling items



PRODUCTS SHOWN ACTUAL SIZE

TAPERED CAPS

Radiopaque and compatible with steam sterilization. For single use.

	REF	COLOR	Ø DIA	QTY
☆ A	31101	white transparent	2.0 - 3.0 x 38mm	70 pcs
B	31201	white	2.0 - 3.0 x 51mm	70 pcs
☆ C	31102	blue transparent	3.0 - 4.0 x 38mm	70 pcs
D	31202	blue transparent	3.0 - 4.0 x 51mm	70 pcs
E	31103	green transparent	4.0 - 5.0 x 51mm	70 pcs
F	31203	green transparent	4.0 - 5.0 x 38mm	70 pcs
G	31104	red transparent	5.0 - 6.0 x 51mm	70 pcs
H	31204	red transparent	5.0 - 6.0 x 38mm	70 pcs
I	31105	yellow transparent	6.0 - 7.0 x 38mm	70 pcs
J	31205	yellow transparent	6.0 - 7.0 x 51mm	70 pcs
K	31106	purple transparent	7.0 - 8.0 x 38mm	70 pcs
L	31206	purple transparent	7.0 - 8.0 x 51mm	70 pcs
M	31107	orange transparent	8.0 - 9.0 x 38mm	70 pcs
N	31207	orange transparent	8.0 - 9.0 x 51mm	70 pcs
O	31108	brown transparent	9.0 - 10 x 38mm	70 pcs
P	31208	brown transparent	9.0 - 10 x 51mm	70 pcs
Q	31109	black transparent	10 - 11 x 38mm	70 pcs
R	31209	black transparent	10 - 11 x 51mm	70 pcs
S	31303	green transparent	13 - 14 x 38mm	50 pcs
T	31403	green transparent	15 - 16 x 51mm	50 pcs
U	31305	yellow transparent	15 - 16 x 38mm	50 pcs
V	31405	yellow transparent	19 - 20 x 51mm	50 pcs
W	31307	orange transparent	15 - 16 x 38mm	50 pcs
X	31407	orange transparent	19 - 20 x 51mm	50 pcs

☆ Top selling items



PRODUCTS SHOWN ACTUAL SIZE

FLAT CAPS

Radiopaque and compatible with steam sterilization. For single use.

	REF	DESCRIPTION	COLOR	Ø DIA	QTY
☆ A	30407	for osteotomes	clear	2.0 x 9.0 x 25mm	100 pcs
B	30807	for osteotomes	orange transparent	2.0 x 9.0 x 25mm	100 pcs
C	31507	for osteotomes	orange	2.0 x 9.0 x 25mm	100 pcs
☆ D	30408	for osteotomes	clear	2.0 x 16 x 25mm	100 pcs
E	30808	for osteotomes	brown transparent	2.0 x 16 x 25mm	100 pcs
F	31508	for osteotomes	brown	2.0 x 16 x 25mm	100 pcs
☆ G	30409	for osteotomes	clear	3.0 x 25 x 25mm	100 pcs
H	30809	for osteotomes	black transparent	3.0 x 25 x 25mm	100 pcs
I	31509	for osteotomes	black	3.0 x 25 x 25mm	100 pcs
☆ J	30411	for osteotomes	clear	6.4 x 19mm	50 pcs
K	30811	for osteotomes	white transparent	6.4 x 19mm	100 pcs
L	30812	for osteotomes	blue transparent	9.5 x 25mm	100 pcs
M	30813	for osteotomes	green	13 x 25mm	100 pcs
N	30815	for osteotomes	yellow	18 x 25mm	100 pcs
O	30415	for osteotomes	clear	19 x 25mm	50 pcs
P	30409	for osteotomes	clear	25 x 25mm	50 pcs
Q	30417	for osteotomes	clear	32 x 32mm	50 pcs
R	30817	for osteotomes	purple transparent	32 x 32mm	100 pcs
S	30418	for osteotomes	clear	38 x 38mm	50 pcs
T	30818	for osteotomes	red	38 x 38mm	100 pcs
U	30419	for osteotomes	clear	45 x 32mm	100 pcs
V	30819	for osteotomes	orange	45 x 32mm	100 pcs

☆ Top selling items



PRODUCTS SHOWN ACTUAL SIZE

SKIN HOOK CAPS

Protect delicate single and double prong skin hook instrumentation. Radiopaque. Compatible with sterilization. Single use.

REF	DESCRIPTION	COLOR	Ø DIA	QTY
A 31003	double prong. 2mm apart	clear	4.0 x 19mm	100 pcs
B 31001	single prong. 2.53.5mm hook	clear	4.8 x 19mm	100 pcs
C 31002	single prong. 5mm hook	clear	6.3 x 19mm	100 pcs
D 31004	double prong. 5mm apart	clear	6.7 x 19mm	100 pcs
E 31005	double prong. 7mm apart	clear	9.1 x 19mm	100 pcs

SELF-RETAINING AND GELPI RETRACTOR CAPS

Protect a variety of self-retaining retractors including Weitlaner, D'Errico-Adson, Beckman, and Adson retractors. Radiopaque. Compatible with steam sterilization. Single use.

REF	DESCRIPTION	COLOR	Ø DIA	QTY
☆ F 30907	selfretaining retractor cap	blue transparent	22 x 20mm	100 pcs
☆ G 30908	selfretaining retractor cap	purple	25 x 25mm	100 pcs
☆ H 30909	selfretaining retractor cap	orange	25 x 32mm	100 pcs
I 30910	selfretaining retractor cap	green	45 x 13mm	100 pcs
J 30903	For Gelpi retractors	clear	51mm	100 pcs
K 30904	For Gelpi retractors	clear	76mm	100 pcs
L 30905	For Gelpi retractors	clear	38mm	100 pcs

☆ Top selling items



PRODUCTS SHOWN ACTUAL SIZE

TWIN CAPS

Protect instruments while keeping them in an open position for cleaning. Perfect for scissors and forceps. Radiopaque. Single use. Compatible with steam sterilization.


	REF	COLOR	Ø DIA	QTY
★ A	30601	white transparent	1.6 x 19mm	50 pcs
★ B	30602	blue transparent	2.0 x 19mm	50 pcs
★ C	30603	green transparent	2.8 x 19mm	50 pcs
★ D	30607	orange transparent	2.0 x 25mm	50 pcs
★ E	30604	red transparent	3.2 x 25mm	50 pcs
★ F	30605	yellow transparent	5.0 x 25mm	50 pcs
	30650	clear	-	1 set

AIRFLOW CAPS™

An alternative to perforated caps, grips on the inside of the cap allow more airflow during sterilization process and assist in easy removal. Radiopaque. Compatible with steam sterilization. Single use.

	REF	COLOR	Ø DIA	QTY
★ G	31608	clear	2.0 x 16 x 25mm	100 pcs
★ H	31605	clear	4.8 x 25mm	100 pcs
★ I	31606	clear	10 x 20mm	100 pcs
★ J	31607	clear	9.5 x 25mm	100 pcs
★ K	31609	clear	25 x 25mm	100 pcs

★ Top selling items




INSTRUMENT ORGANIZER

Designed to help protect and organize your delicate surgical instruments within a wrapped set or container system. Single use. Compatible with steam sterilization.

REF	DESCRIPTION	LENGTH	WIDTH	QTY
A	CSDH100 2 rows of 6 pockets	444.5mm	292.1mm	100 pcs
B	CSDH102 2 rows of 6 pockets with slots in pocket	444.5mm	292.1mm	100 pcs
C	CSDH200 2 rows of 6 notched pockets with tape closure strips	470mm	260mm	100 pcs
D	CSDH300 2 rows of 6 pockets	460mm	410mm	100 pcs
E	CSDH505 single pocket	140mm	133mm	100 pcs
F	CSDH512 single pocket	127mm	305mm	100 pcs

☆
☆
☆
☆



☆ Top selling items

TAGS

Our tags allow you to indicate special handling instructions, secure trays and carts from tampering and quickly identify priority instruments for repair and sterilization. Tags are available in a wide range of styles, colors and locking configurations. Printed tags with special handling instructions, as well as uniquely numbered tags, are available as a standard product.

Write-On Tags allow you to write specific handling instructions on each tag when more detailed communication is necessary. These tags attach to instruments and trays with ties, and are larger and more visible than other tags, ensuring that detailed instructions will be seen easily.

Metal Tray Tags clip onto the edge of a tray and provide ample space for identification, tracking and barcoding needs. Made with stainless steel for durability.

[Plastic Tray Tags](#)

[Clamp Tag Holders](#)

[Metal Tray Tags](#)

[Tags Accessories](#)

[One-Loop Tags](#)

[Self-Tie Tags](#)

[Write-On Tag](#)






FIRST TAG SHOWN ACTUAL SIZE


TRAY TAGS FOR COLOR CODING

Durable, plastic tags designed for use on instrument trays. Tray tags provide a flat surface for barcode placement. Great for color coding. Compatible with steam and H₂O₂ sterilization. Designed to be attached to a tray using a fastener such as a ring or clip.

	REF	COLOR	FORMAT	QTY
A	01900	red	90 x 33mm	100 pcs
B	01901	green	90 x 33mm	100 pcs
C	01902	light blue	90 x 33mm	100 pcs
D	01904	white	90 x 33mm	100 pcs
E	01905	orange	90 x 33mm	100 pcs
F	01912	black	90 x 33mm	100 pcs
G	01956	brown	90 x 33mm	100 pcs
H	01962	dark blue	90 x 33mm	100 pcs
I	01976	pink	90 x 33mm	100 pcs
J	01977	purple	90 x 33mm	100 pcs



FIRST TAG SHOWN ACTUAL SIZE



TAGS FOR BASKET IDENTIFICATION

Durable, plastic tags designed for use on instrument trays. Tray tags provide a flat surface for barcode placement. Great for color coding. Compatible with steam and H₂O₂ sterilization. The tags are designed to be attached to a tray using a fastener such as a ring or clip.

	REF	COLOR	FORMAT	QTY
☆ A	01948	white	60 x 40mm	100 pcs
☆ B	01949	red	60 x 40mm	100 pcs
☆ C	01950	green	60 x 40mm	100 pcs
D	01969	dark blue	60 x 40mm	100 pcs
E	01951	light blue	60 x 40mm	100 pcs
F	01952	orange	60 x 40mm	100 pcs
G	01953	black	60 x 40mm	100 pcs
☆ H	01954	yellow	60 x 40mm	100 pcs
I	01990	pink	60 x 40mm	100 pcs
J	01991	purple	60 x 40mm	100 pcs
K	02278	brown	60 x 40mm	100 pcs
L	02279	ochre yellow	60 x 40mm	100 pcs

☆ Top selling items



PLASTIC TRAY TAGS

Tray tags simplify the identification of sets or instruments and provide a useful means of labeling items for tracking and tracing. All tags are produced to withstand the rigorous decontamination cycle. Available in 18 colors. Tray tags can be attached to trays using either a split ring or a cable tie.

	REF	COLOR	FORMAT	QTY
A	TAG020	black	90 x 38mm	100 pcs
B	TAG018	gray	90 x 38mm	100 pcs
C	TAG026	royal blue	90 x 38mm	100 pcs
D	TAG023	turquoise	90 x 38mm	100 pcs
☆ E	TAG017	green	90 x 38mm	100 pcs
F	TAG021	lime	90 x 38mm	100 pcs
☆ G	TAG013	yellow	90 x 38mm	100 pcs
H	TAG014	orange	90 x 38mm	100 pcs
☆ I	TAG015	red	90 x 38mm	100 pcs
J	TAG012	brown	90 x 38mm	100 pcs
☆	TAG011	white	90 x 38mm	100 pcs
☆	TAG016	blue	90 x 38mm	100 pcs
	TAG019	purple	90 x 38mm	100 pcs
	TAG022	pink	90 x 38mm	100 pcs
	TAG024	mustard	90 x 38mm	100 pcs
	TAG025	terracotta	90 x 38mm	100 pcs

☆ Top selling items



FIRST TAG SHOWN ACTUAL SIZE

PTL TIE-ON TAGS

PTL tie-on tags are designed for multiple use. They can be attached to the tray by using REF 01589. These are compatible with P-touch and laser inscribed including text and barcode. Heat-resistant up to 134 °C.

	REF	COLOR	FORMAT	QTY
A	02253	light blue	90 x 40mm	100 pcs
B	02254	white	90 x 40mm	100 pcs

MESH TRAY TAGS

Mesh tray tags are not suitable for P-touch, and laser marking is not possible. Heat-resistant up to 134 °C.

	REF	DESCRIPTION	FORMAT	QTY
C	01421	metallic	90 x 50mm	150 pcs
	01422	metallic	90 x 33mm	200 pcs






PLASTIC TRAY TAGS

Durable, plastic tags designed for use on instrument trays. Tray tags provide a flat surface for barcode placement. Great for color coding. Compatible with steam and H₂O₂ sterilization. Tags include a built-in clip on the back.

	REF	COLOR	FORMAT	QTY
A	01563	yellow	30 x 30mm	25 pcs
B	01913	white	30 x 30mm	25 pcs
C	01914	red	30 x 30mm	25 pcs
D	01915	green	30 x 30mm	25 pcs
E	01916	light blue	30 x 30mm	25 pcs
F	01917	orange	30 x 30mm	25 pcs
G	01919	ochre yellow	30 x 30mm	25 pcs
H	01961	brown	30 x 30mm	25 pcs
I	01981	purple	30 x 30mm	25 pcs
J	01561	yellow	55 x 30mm	25 pcs
K	01920	white	55 x 30mm	25 pcs
L	01921	red	55 x 30mm	25 pcs
M	01922	green	55 x 30mm	25 pcs
N	01923	light blue	55 x 30mm	25 pcs
O	01924	orange	55 x 30mm	25 pcs
P	01925	black	55 x 30mm	25 pcs
Q	01926	ochre yellow	55 x 30mm	25 pcs
R	01965	dark blue	55 x 30mm	25 pcs
S	01982	pink	55 x 30mm	25 pcs
T	01983	purple	55 x 30mm	25 pcs

FIRST TAG SHOWN ACTUAL SIZE







FIRST TAG SHOWN ACTUAL SIZE

PLASTIC TRAY TAGS

Durable, plastic tags designed for use on instrument trays. Tray tags provide a flat surface for barcode placement. Great for color coding. Compatible with steam and H₂O₂ sterilization. Tags include a built-in clip on the back.

	REF	COLOR	FORMAT	QTY
A	01927	white	70 x 30mm	25 pcs
B	01928	red	70 x 30mm	25 pcs
C	01929	green	70 x 30mm	25 pcs
D	01930	light blue	70 x 30mm	25 pcs
E	01931	orange	70 x 30mm	25 pcs
F	01932	black	70 x 30mm	25 pcs
G	01933	ochre yellow	70 x 30mm	25 pcs
H	01966	dark blue	70 x 30mm	25 pcs
I	01985	purple	70 x 30mm	25 pcs
J	01564	yellow	90 x 30mm	25 pcs
K	01934	white	90 x 30mm	25 pcs
L	01935	red	90 x 30mm	25 pcs
M	01936	green	90 x 30mm	25 pcs
N	01937	light blue	90 x 30mm	25 pcs
O	01938	orange	90 x 30mm	25 pcs
P	01939	black	90 x 30mm	25 pcs
Q	01940	ochre yellow	90 x 30mm	25 pcs
R	01958	brown	90 x 30mm	25 pcs
S	01967	dark blue	90 x 30mm	25 pcs



FIRST TAG SHOWN ACTUAL SIZE

PTL INSERT LABELS

PTL insert labels are reusable. They are suitable for REF 01548, 01497, 01498 and 01995 as well as container holders. These are P-touch compatible and laser labeling including text and barcode is possible. Compatible with steam sterilization and heat-resistant up to 134 °C.

	REF	COLOR	FORMAT	QTY
A	01566	white	55 x 35mm	25 pcs
B	01567	red	55 x 35mm	25 pcs
C	01568	green	55 x 35mm	25 pcs
D	01569	light blue	55 x 35mm	25 pcs
E	01570	yellow	55 x 35mm	25 pcs
F	01571	orange	55 x 35mm	25 pcs
G	01573	dark blue	55 x 35mm	25 pcs
H	01575	pink	55 x 35mm	25 pcs
I	01576	purple	55 x 35mm	25 pcs
J	01550	white	56 x 18mm	25 pcs
K	01551	red	56 x 18mm	25 pcs
L	01552	green	56 x 18mm	25 pcs
M	01553	light blue	56 x 18mm	25 pcs
N	01554	yellow	56 x 18mm	25 pcs
O	01555	orange	56 x 18mm	25 pcs
P	01556	black	56 x 18mm	25 pcs
Q	01970	dark blue	56 x 18mm	25 pcs
R	01992	pink	56 x 18mm	25 pcs
S	01993	purple	56 x 18mm	25 pcs
T	01547	gray	56 x 18mm	25 pcs
U	01971	coral red	56 x 18mm	25 pcs
V	02215	white	70 x 35mm	25 pcs
W	02216	red	70 x 35mm	25 pcs
X	02217	green	70 x 35mm	25 pcs
Y	02219	dark blue	70 x 35mm	25 pcs
Z	02220	yellow	70 x 35mm	25 pcs
AA	01598	black	70 x 35mm	25 pcs
BB	02221	white	90 x 35mm	25 pcs
CC	02222	red	90 x 35mm	25 pcs
DD	02223	green	90 x 35mm	25 pcs
EE	02224	light blue	90 x 35mm	25 pcs
FF	02225	dark blue	90 x 35mm	25 pcs
GG	02226	yellow	90 x 35mm	25 pcs
HH	01599	black	90 x 35mm	25 pcs



CLAMP HOLDER MADE OF STAINLESS STEEL

Durable and reusable holders made of stainless steel and heat-resistant.


	REF	DESCRIPTION	FORMAT	QTY
A	01548	tag holder with side flap, suitable for PTL insertion labels	60 x 40mm	25 pcs
B	01497	tag holder, suitable for PTL insert and aluminum tags	50 x 20mm	25 pcs
C	01498	tag holder, suitable for PTL tag	60 x 40mm	25 pcs

ALUMINUM LABELS

Suitable for narrow inserts on containers and laser marking.

	REF	COLOR	FORMAT	QTY
D	28500	silver	18 x 50mm	25 pcs
E	28501	blue	18 x 50mm	25 pcs
F	28502	red	18 x 50mm	25 pcs






TAGS SHOWN ACTUAL SIZE

TAGS MADE OF STAINLESS STEEL FOR SIEVE TRAY LABELING

Durable tags made of stainless steel. Reusable and heat-resistant.

	REF	DESCRIPTION	FORMAT	QTY
A	01444	label	90 x 50mm	50 pcs
B	01445	label	90 x 33mm	50 pcs
C	01452	tag with hole (4 x 9mm)	Ø 25mm	50 pcs
D	01453	tag with hole (4 x 9mm)	Ø 50mm	50 pcs
E	01588	tie-on tag with hole (6 x 10mm)	65 x 24mm	50 pcs
☆ F	01449	clip label	90 x 30mm	50 pcs
G	01454	clip label	90 x 50mm	50 pcs
H	01559	clip label	55 x 30mm	50 pcs
☆ I	01493	clip label	30 x 30mm	50 pcs
J	01590	clip label, suitable for double frame	75 x 35mm	50 pcs
K	01520	clip label, suitable for double frame	42 x 36mm	50 pcs



☆ Top selling items

A



STOPPERS FOR TRAY ACCESSORIES

Stoppers designed to secure tray accessories such as labels, tags, corner protectors, among others.

A

REF	MATERIAL	COLOR	QTY
09097	silicone	uncolored	100 pcs

ACCESSORIES TO SECURE TRAY ACCESSORIES



B

C

REF	DESCRIPTION	FORMAT	QTY
TG100	clipper for tie-on tag; with holes	40mm	5 pcs
TG200	key ring for affixing sieve tray tags	Ø 20mm	100 pcs


B



C



☆ Top selling items




PLASTIC TRAY LABELS BLANK

Designed to loop around tray handles to identify through color coding or providing a flat surface for barcode placement. Compatible with steam and EO sterilization. Single use.

	REF	COLOR	FORMAT	QTY
☆ A	09634	white	160 x 35mm	250 pcs
☆ B	09637	red	160 x 35mm	250 pcs
☆ C	09636	green	160 x 35mm	250 pcs
☆ D	09638	yellow	160 x 35mm	250 pcs
E	09676	brown	160 x 35mm	250 pcs
F	09635	blue	160 x 35mm	250 pcs
☆ G	09268	white	206 x 110mm	250 pcs
	09751	white	240 x 50mm	250 pcs

☆ Top selling items



PLASTIC TRAY LABELS BLANK

For single use, heat-resistant up to 134 °C, suitable for washer-disinfectors.

	REF	COLOR	FORMAT	QTY
A	09264	blue	236 x 50mm	250 pcs
	09566	blue one-sided	240 x 50mm	250 pcs
B	09267	yellow	236 x 50mm	250 pcs
	09569	yellow one-sided	240 x 50mm	250 pcs
C	09265	green	236 x 50mm	250 pcs
	09567	green one-sided	240 x 50mm	250 pcs
D	09269	orange	236 x 50mm	250 pcs
E	09266	red	236 x 50mm	250 pcs
	09568	red one-sided	240 x 50mm	250 pcs
F	09244	white	236 x 50mm	250 pcs
G	08338	pink one-sided	240 x 50mm	250 pcs





SELF-TIE TAGS


Self-Tie Tags are designed for the identification of trays that are to be steam sterilized. Plain tag with dumbbell self-tie. Single use.

REF	DESCRIPTION	FORMAT	QTY
09691	reel 5,000 pcs	160 x 19mm	1 roll

POLYPROPYLENE TRAY TAGS

Compatible with steam sterilization and heat-resistant up to 134 °C, suitable for washer-disinfectors.

REF	COLOR	LETTERING FIELD SIZE	LENGTH	QTY
09679	white	50 x 30mm	200mm	100 pcs
09681	green	50 x 30mm	200mm	100 pcs
09683	yellow	50 x 30mm	200mm	100 pcs



A



B



C



D



E



F



G



DOUBLE-SIDED

WRITE-ON TAGS

Durable and economical. Identify and record specific instructions with these versatile tags. Tamper-resistant. Compatible with steam and H₂O₂ sterilization. Customizable. 76 x 127mm. Each package includes 100 cable ties and a lab marker. Single use.

	REF	DESCRIPTION	COLOR	QTY
A	WT0009	BLANK	black	100 pcs
B	WT0102	VENDOR	blue	100 pcs
C	WT0203	LOANER	green	100 pcs
D	WT0304	PRIORITY	red	100 pcs
E	WT0405	RERUN	yellow	100 pcs
F	WT0507	BROKEN	orange	100 pcs
G	WT1911	MISSING	pink	100 pcs
H	RT100	BROKEN DO NOT USE	white	100 pcs
I	RT200	BLANK	white	100 pcs

H




I





A

DOUBLE-SIDED



B

DOUBLE-SIDED

WRITE-ON TAGS

Durable and economical. Identify and record specific instructions with these versatile tags. Tamper-resistant. Compatible with steam and H₂O₂ sterilization. Customizable. 76 x 127mm. Each package includes 100 cable ties and a lab marker. Single use.

	REF	DESCRIPTION	COLOR	QTY
A	WT2614	LOANER	light green	100 pcs
B	WT3217	QUICK TURNOVER	blue	100 pcs

TRAYS AND ACCESSORIES

Strict requirements have to be met when it comes to the reprocessing of various instruments and medical devices. We offer a wide range of these for optimizing individual procedures and keeping instrument turnaround times to a minimum.

In order to meet the needed requirements, we offer mesh and sterilization trays in a variety of designs, sizes and made of premium material for versatile use (e.g., ophthalmology).

In our product portfolio you will also find comprehensive mesh tray accessories for multiple applications and to protect your sensitive instruments.

Mesh Trays

Special Mesh and Plastic Trays

Instrument Clips

Silicone Mats

Grid Inserts and Mats

Corner Protectors



A



We supply a wide range of stainless steel sieve trays which satisfy the highest of standards in order to guarantee optimal and efficient processes. Our mesh trays made in Germany (Kögel mesh baskets) with a long shelf life are available in DIN and ISO formats in the most widely used sizes, with different mesh grid size and designs.

MESH TRAY WITH DROP HANDLES

Mesh trays without a lid. Made of high quality stainless steel.

REF	FORMAT	MESH SIZE	QTY
08909	250 x 240 x 50mm	6.0 x 6.0mm	1 pc
08951	280 x 160 x 30mm	5.0 x 5.0mm	1 pc
08908	400 x 250 x 70mm	6.0 x 6.0mm	1 pc

A

B

B








SMALL PARTS TRAYS WITH LID

Mesh trays with lid. Made of stainless steel.

	REF	FORMAT	MESH SIZE	QTY
☆ A	08910	40 x 40 x 25mm	1.0 x 1.0mm	1 pc
B	08896	40 x 40 x 25mm	3.0 x 3.0mm	1 pc
C	08897	80 x 40 x 20mm	3.0 x 3.0mm	1 pc
☆ C	08935	80 x 40 x 25mm	1.0 x 1.0mm	1 pc
D	08435	80 x 40 x 30mm	1.6 x 1.6mm	1 pc
E	08900	80 x 40 x 40mm	3.0 x 3.0mm	1 pc
F	08901	80 x 80 x 40mm	3.0 x 3.0mm	1 pc
	08902	90 x 90 x 50mm	3.0 x 3.0mm	1 pc
	08903	140 x 110 x 20mm	3.0 x 3.0mm	1 pc
	08911	80 x 80 x 40mm	1.0 x 1.0mm	1 pc
☆ G	08931	150 x 110 x 20mm	1.0 x 1.0mm	1 pc
	08912	160 x 40 x 20mm	1.0 x 1.0mm	1 pc
☆ H	08898	160 x 40 x 20mm	3.0 x 3.0mm	1 pc
☆ I	08904	180 x 130 x 20mm	3.0 x 3.0mm	1 pc

☆ Top selling items



SMALL PARTS TRAYS

Mesh trays with lid. Made of stainless steel.

	REF	FORMAT	MESH SIZE	QTY
☆ A	08936	200 x 40 x 40mm	1.0 x 1.0mm	1 pc
B	08899	200 x 40 x 40mm	3.0 x 3.0mm	1 pc
C	08948	200 x 60 x 40mm	3.0 x 3.0mm	1 pc
D	08772	210 x 80 x 40mm	5.0 x 5.0mm	1 pc
E	08905	210 x 170 x 20mm	3.0 x 3.0mm	1 pc
☆ F	08957	230 x 150x 50mm	1.0 x 1.0mm	1 pc
	08906	250 x 200 x 20mm	3.0 x 3.0mm	1 pc
	08875	300 x 150 x 70mm	1.0 x 1.0mm	1 pc
G	08950	300 x 150 x 70mm	4.0 x 4.0mm	1 pc

☆ Top selling items



KÖGEL SYSTEM TRAYS AND SETS

Mesh trays with lid. Made of chromium-nickel steel 18/10 (1.4301), electrolytically polished, burr-free, with double peripheral safety frame (DBGM).

REF	DESCRIPTION	DIN STE STANDARD	FORMAT	MESH SIZE	QTY
A 08400	1 x system small parts tray with clamps lid for variable partitioning: 1 x lengthwise divider 1/1 1 x crosswise divider 1/1 1 x crosswise divider 2/3 1 x crosswise divider 1/3 8 x T-clamps	1/1	480 x 250 x 57mm	2.5 x 2.5 x 0.7	1 set
B 08469	1 x system small parts tray with clamps lid	1/3	466 x 74 x 28mm	1.0 x 1.0 x 0.5	1 set



KÖGEL SYSTEM

Mesh trays without lid.lid. Made of chromium-nickel steel.

REF	DIN STE STANDARD	FORMAT	MESH SIZE	QTY
A 08415	1/2	240 x 250 x 50mm	4.0 x 4.0mm	1 pc
B 08416	1/1	480 x 250 x 50mm	4.0 x 4.0mm	1 pc
08809	XL	560 x 250 x 50mm	4.0 x 4.0mm	1 pc
08810	XL	560 x 250 x 70mm	5.0 x 5.0mm	1 pc

KÖGEL SYSTEM SMALL PARTS

Wire container with lid. Made of chromium-nickel steel electrolytically polished, burr-free, with a safe double frame (DBGM).

REF	DESCRIPTION	DIN STE STANDARD	FORMAT	MESH SIZE	QTY
08463	set with lid and clamps	1/4	236 x 121 x 28mm	1.0 x 1.0 x 0.5mm	1 set
08404	set with lid and clamps	1/4	236 x 121 x 46mm	2.5 x 2.5 x 0.7mm	1 set
08454	set with lid and clamps	1/6	236 x 121 x 28mm	1.0 x 1.0 x 0.5mm	1 set
08457	without lid and clamps	1/6	236 x 73.5 x 28mm	1.0 x 1.0 x 0.5mm	1 pc



A



B




C



KÖGEL SYSTEM SMALL PARTS

Mesh trays with lid. Made of chromium-nickel steel electrolytically polished, burr-free, and double peripheral safety frame (DBGM).

REF	DESCRIPTION	DIN STE STANDARD	FORMAT	MESH SIZE	QTY
A 08406	set with lid and clamps	1/8	121 x 121 x 29mm	2.5 x 2.5 x 0.7mm	1 set
B 08408	set with lid and clamps	1/16	121 x 73 x 29mm	1.0 x 1.0 x 0.5mm	1 set
C 08442	set with lid and clamps	1/32	63 x 63 x 23.4mm	1.0 x 1.0 x 0.5mm	1 set



MESH TRAYS FOR SPECIAL APPLICATIONS


REF	DESCRIPTION	LID	FORMAT	QTY
08774	Mesh basket with lid for Da Vinci arms, mesh size: 5.5 x 5.5 mm, made of stainless steel, suitable for: Silicone holder REF A327593, suitable for Da Vinci systems 470006 and 470179	yes	635 x 95 x 60mm	1 pc
645DVXL	Mesh tray for Da Vinci XI system	no	640 x 340 x 100mm	1 pc
79DVXL1	sterilization tray for rigid-flexible endoscope	yes	670 x 80 x 52mm	1 pc

STERILIZATION TRAYS AND CASSETTES

Trays and cassettes exclusively for sterilization and instrument storage. Made of plastic.

REF	DESCRIPTION	FORMAT	HEAT RESISTANCE	QTY
08892	with lid and silicone studded mat	185 x 75 x 20mm	up to 134 °C	1 pc
08890	with lid and silicone studded mat	190 x 100 x 20mm	up to 134 °C	1 pc





INSTRUMENT CLIPS


These clips can be used to bundle ring instruments so that instrument groups are stored together tidily in the mesh tray. Made of stainless steel.


	REF	TOTAL LENGTH	QTY
A	20403	75mm	1 pc
B	20404	100mm	1 pc
C	20405	135mm	1 pc

SILICONE INSERTS

With pins for fastening in the tray, ideally suited for a mesh size of 5mm x 5mm to 6mm x 6mm.

	REF	DESCRIPTION	FORMAT	Ø DIA	QTY
D	08949	round opening (5-slots)	125 x 25mm	4mm	6 pcs





SILICONE NUB MATS


Silicone mats designed to provide support to instruments during high level disinfection. Suitable for washer-disinfectors, flexible.

REF	COLOR	FORMAT	HEAT RESITSTANCE	QTY
08698	pink	245 x 1000mm	up to 230 °C	1 pc
08321	blue	500 x 320mm	up to 40 °C	1 pc
08609	blue	470 x 250 x 3.5mm	up to 230 °C	1 pc


☆ **A**
B

☆ Top selling items


A



B




C



D






E



GRID INSERTS AND MATS

Silicone grids and mats designed to provide support to instruments during high level disinfection. Suitable for washer-disinfectors, honeycomb structure and flexible.

	REF	DESCRIPTION	COLOR	FORMAT	HEAT RESISTANCE	QTY
A	08602	honeycomb structure	dark blue	238 x 248mm	up to 230 °C	1 pc
B	08608	honeycomb structure	light green	270 x 120mm	up to 230 °C	1 pc
C	08601	honeycomb structure	dark blue	480 x 250mm	up to 230 °C	1 pc
D	08110	silicone grid insert, thickness 6mm	white	600 x 324mm	up to 134 °C	1 pc
E	08605	silicone grid insert, thickness 2mm	blue	600 x 400mm	-40° C up to +230°	1 pc
	08600	silicone grid insert, thickness 2mm	blue	540 x 320mm	-60° C to +300 °C	1 pc



CORNER PROTECTORS

Corner protectors are designed to protect the soft packaging from the sharp edges of the sterilization tray.

	REF	DESCRIPTION	COLOR	FORMAT	QTY
☆ A	CS9113	silicone corner protector	lime green	71.12 x 71.12 x 30.48mm	80 pcs
☆ B	CS9125C	foam corner protector for single use	blue	152.4 x 152.4 x 31.75mm	80 pcs
☆ C	CC002	crepe corner protector for single use	white	200 x 150mm	100 pcs
☆ D	CC001	cardboard corner protector for single use	white	200 x 150mm	200 pcs
	42007	cardboard corner protector for single use	white	100 x 100 x 50mm	600 pcs
☆	42010	cardboard corner protector for single use	white	100 x 100 x 50mm	600 pcs

☆ Top selling items

CONTAINER ACCESSORIES

Containers are used worldwide to provide a protective barrier during the sterilization, handling and storage of surgical instruments. Regardless of which container model you are using, we have just the right accessories.

Container filters are available in various shapes and sizes for a variety of container systems in accordance with the applicable standards.

Colored plastic seals in various shapes that can be used for color coding and as a sustainable container closure option, paper seals are perfect for easy handling.

Container Filters

Plastic Locking Seals

Numbered Locking Seals

Paper Locking Seals

Insert Labels

Tag Labels




CONTAINER FILTERS


We offer a variety of filters for containers compliant with DIN EN ISO 11607 standard. These paper filters are available with or without process indicator and come in different sizes and shapes. Single use.

CONTAINER FILTERS FOR AESCULAP CONTAINERS


REF	DESCRIPTION	FORMAT	QTY
A	01152	round filter with steam indicator, central hole punching Ø 6mm	Ø 190mm
B	01141	round filter without indicator, central hole punching Ø 6mm, white	Ø 190mm
C	01148	small set filter with steam indicator, rectangular, white	118 x 235mm
D	01186	small set filter without indicator, rectangular, white	118 x 235mm
E	01317	rectangular filter with steam indicator	170 x 230mm
F	01153	square filter with steam indicator, white	230 x 230mm
G	01142	square filter without indicator, white	230 x 230mm
H	01149	square filter without indicator, white	230 x 230mm
I	01330	square filter with steam indicator	104 x 104mm
J	01158	rectangular filter with steam indicator	110 x 240mm
K	01143	rectangular filter with steam indicator	180 x 270mm
L	01321	rectangular filter with steam indicator	180 x 255mm

Top selling items





A



B



C

CONTAINER FILTERS

Rectangular, white. Compliant with DIN EN ISO 11607 standard. Single use.

	REF	DESCRIPTION	FORMAT	QTY
A	01144	rectangular, without indicator	240 x 240mm	250 pcs
	01198	rectangular, without indicator	230 x 540mm	250 pcs
B	01304	rectangular, without indicator	210 x 490mm	100 pcs

NON-WOVEN CONTAINER FILTERS

Round and square format container filters, 67gsm according to DIN EN ISO 11607. Single use.

	REF	DESCRIPTION	FORMAT	QTY
☆ C	01145	square, green, with steam indicator,	225 x 225mm	500 pcs
☆ D	01378	round, blue, for H ₂ O ₂ sterilization, with steam indicator	190 x 6.0mm	500 pcs
	01365	round, blue, with steam indicator	190 x 6.0mm	500 pcs
E	01368	round, green, without indicator, no filter imprint	Ø 92mm	1000 pcs




D



E


☆ Top selling items



A



B



C




D



E



F



G

CONTAINER FILTERS

For single use.

REF	DESCRIPTION	FORMAT	QTY
01314	rectangular filter, with indicator	145 x 205mm	100 pcs
01315	square filter, with indicator	240 x 240mm	250 pcs
01327	rectangular filter, with indicator	240 x 540mm	100 pcs

A

B

C

CONTAINER FILTERS FOR WAGNER CONTAINERS

Round and rectangular, white, with indicator. Single use.

REF	DESCRIPTION	FORMAT	QTY
01355	round, central hole punching Ø 10mm	Ø 190mm	500 pcs
01369	round, central hole punching Ø 10mm	Ø 75mm	1000 pcs
01363	rectangular	245 x 400mm	100 pcs
01323	rectangular	96 x 212mm	250 pcs


D

E

F

G

☆ Top selling items



PLASTIC LOCKING SEALS





We offer a wide range of different designs to suit all type of containers. Our seals are available both with and without process indicators. They are used as a tamper free lock and have a predetermined breaking point for easy opening. We have options for all common container makes and a wide range of colors. We also offer numbered security seals (REF 02639) for the sealed return of loan instruments, implants, and containers. Single use.

LOCKING SEAL - ANGULAR BRIDGE

Bridge length 19mm.

	REF	DESCRIPTION	COLOR	QTY
☆ A	01020	without indicator	dark blue	1000 pcs
☆ B	01021	without indicator	green	1000 pcs
☆ C	01022	without indicator	white	1000 pcs
☆ D	01023	without indicator	red	1000 pcs
☆ E	01024	without indicator	yellow	1000 pcs
☆ F	01234	without indicator	orange	1000 pcs
☆ G	01286	without indicator	black	1000 pcs
☆ H	02623	without indicator	light blue	1000 pcs
☆ I	01228	with steam indicator	dark blue	1000 pcs
☆ J	01229	with steam indicator	green	1000 pcs
☆ K	01230	with steam indicator	white	1000 pcs
☆ L	01231	with steam indicator	red	1000 pcs
☆ M	01232	with steam indicator	yellow	1000 pcs
☆ N	01235	with steam indicator	orange	1000 pcs
☆ O	02670	with steam indicator	light blue	1000 pcs

☆ Top selling items



UNIVERSAL LOCKING SEALS UNILOCK - ANGULAR BRIDGE




Bridge length 19mm. Single use.

	REF	DESCRIPTION	COLOR	QTY
A	01240	without indicator	blue	1000 pcs
B	01241	without indicator	green	1000 pcs
C	01242	without indicator	white	1000 pcs
D	01246	with steam indicator	green	1000 pcs
E	01248	with steam indicator	red	1000 pcs
F	01249	with steam indicator	yellow	1000 pcs

UNIVERSAL LOCKING SEALS UNILOCK XL - SQUARE BRIDGE

Bridge length 27mm. Single use.

	REF	DESCRIPTION	COLOR	QTY
G	01250	without indicator	blue	1000 pcs
H	01251	without indicator	green	1000 pcs
I	01252	without indicator	white	1000 pcs
J	01254	without indicator	yellow	1000 pcs
K	02514	without indicator	violet	1000 pcs
L	02515	without indicator	orange	1000 pcs
M	01255	with steam indicator	blue	1000 pcs
N	01256	with steam indicator	green	1000 pcs
O	01257	with steam indicator	white	1000 pcs
P	01258	with steam indicator	red	1000 pcs
Q	01259	with steam indicator	yellow	1000 pcs



UNIVERSAL LOCKING SEALS - ANGULAR BRIDGE

Suitable for Aesculap and Wagner containers. Bridge length 34mm. Single use.


	REF	DESCRIPTION	COLOR	QTY
A	01025	without indicator	dark blue	1000 pcs
B	01026	without indicator	red	1000 pcs
C	01027	without indicator	green	1000 pcs
D	01028	without indicator	yellow	1000 pcs
E	01029	without indicator	white	1000 pcs
F	02624	without indicator	orange	1000 pcs
G	02625	without indicator	violet	1000 pcs
H	01220	with steam indicator	dark blue	1000 pcs
I	01221	with steam indicator	red	1000 pcs
J	01233	with steam indicator	green	1000 pcs
K	01238	with steam indicator	yellow	1000 pcs
L	01239	with steam indicator	white	1000 pcs

'AW' SEALS - ROUND BAR

Suitable for Aesculap and Wagner containers. Bridge length 19mm. Single use.

	REF	DESCRIPTION	COLOR	QTY
☆ M	02612	without indicator	red	1000 pcs
☆ N	02613	without indicator	blue	1000 pcs
☆ O	02614	without indicator	green	1000 pcs
P	02615	without indicator	yellow	1000 pcs
Q	02661	without indicator	white	1000 pcs
R	02662	without indicator	black	1000 pcs
S	02632	with steam indicator	red	1000 pcs
T	02633	with steam indicator	blue	1000 pcs
U	02634	with steam indicator	green	1000 pcs
V	02635	with steam indicator	yellow	1000 pcs

☆ Top selling items



UNIVERSAL SEALS - ROUND BRIDGE

Length 195mm. Single use.

	REF	DESCRIPTION	COLOR	QTY
A	01204	without indicator	dark blue	1000 pcs
B	02643	without indicator	yellow	1000 pcs
C	01205	with steam indicator	dark blue	1000 pcs

UNIVERSAL SEALS - ANGULAR BRIDGE

Bridge length 80mm. Single use.

	REF	DESCRIPTION	COLOR	QTY
D	01014	without indicator	blue	1000 pcs
E	01216	without indicator	red	1000 pcs
F	01217	without indicator	green	1000 pcs
G	01218	without indicator	yellow	1000 pcs
H	01219	without indicator	black	1000 pcs
☆ I	01013	with steam indicator	blue	1000 pcs
☆ J	01212	with steam indicator	red	1000 pcs
K	01214	with steam indicator	yellow	1000 pcs

☆ Top selling items



A




B

PLASTIC SEALS FOR MARTIN CONTAINERS

REF	DESCRIPTION	COLOR	QTY
A 01011	without indicator	white	1000 pcs
B 01012	with steam indicator	white	1000 pcs

SECURITY SEALS

REF	DESCRIPTION	COLOR	QTY
C 02639	290mm long, usable length 220mm	yellow	400 pcs



C



A

B

C



D





E

LOCKING TAGS

Versatile, durable, and economical. Tamper-resistant. Compatible with steam and H₂O₂ sterilization. Single use. Customizable. 25.4mm length.

	REF	DETAILS	COLOR	QTY
A	LT100N	uniquely numbered	white	100 pcs
	LT100N1000	uniquely numbered	white	1000 pcs
☆ B	LT300N	uniquely numbered	green	100 pcs
	LT300N1000	uniquely numbered	green	1000 pcs
☆ C	LT400N	uniquely numbered	red	100 pcs
	LT400N1000	uniquely numbered	red	1000 pcs
D	LT500N	uniquely numbered	yellow	100 pcs
	LT500N1000	uniquely numbered	yellow	1000 pcs
E	LT900N	uniquely numbered	black	100 pcs
	LT900N1000	uniquely numbered	black	1000 pcs

☆ Top selling items



PAPER SEALS

We offer you an environmentally friendly way of closing containers. The seals can be attached to almost all common containers and follow a simple principle. With our paper seals you ensure safe closure for transport and storage due to thermal closure of the components. Available in different versions with and without documentation fields. In some cases, it is also possible to stick on labels. Single use.

	REF	DESCRIPTION	COLOR	FORMAT	QTY
A	01010	bridge width 3mm	white	100 x 18mm	1000 pcs
B	01003	with indicator, bridge width 3mm	white	105 x 30mm	1000 pcs
C	01004	with indicator, with grip tab, bridge width 4mm	white	117 x 38mm	500 pcs
D	01005	with pull tab, self-adhesive	white	117 x 38mm	500 pcs
E	02638	with indicator, bridge width 2mm	white	105 x 30mm	1000 pcs
F	04120	with indicator, with grip tab, bridge width 4mm	white	132 x 38mm	500 pcs
G	04179	with indicator, self-adhesive	white	145 x 35mm	750 pcs
H	01006	with indicator, with documentation field, self-adhesive, insertion from the right	white	35 x 74mm	1000 pcs
	04114	see 01006, blank	white	35 x 74mm	1000 pcs
I	01008	with indicator, with documentation field, not self-adhesive, insertion from the right	white	35 x 74mm	1000 pcs
J	01009	with indicator, with documentation field, self-adhesive, insertion from the left	white	35 x 74mm	1000 pcs
K	04108	with indicator, with documentation field, not self-adhesive, insertion from the left	white	35 x 74mm	1000 pcs
L	04113	with indicator, with documentation field, self-adhesive, insertion from the right	white	35 x 80mm	1000 pcs
M	61006	closure label, with indicator, self-adhesive, insert from the right, with documentation field	white	74 x 35mm	1000 pcs

☆ Top selling items

REF 01133

ster. no. _____ load. no. _____

ster. date _____

exp. date _____

name _____

steam ↕ brown

A

DAMPF = BRAUN STEAM = BROWN

Steri.-Datum / Date: _____

Verfall / Expiration: _____

Autoclave #: _____ Charge / Lot: _____

Name: _____

REF 01134

B

Date: _____

Exp. Date: _____

Autoclave #: _____ Lot: _____

Name: _____

STEAM = BROWN

REF 01169

C

REF 04067

Steri-Datum / Date: _____

Verfall / Valid until: _____

Name: _____

DAMPF = BRAUN STEAM = BROWN

D

STER. DAT. _____

STER. DATE _____

VERFALL-DAT. _____

EXP. DATE _____

STER. NR. _____

STER. NO. _____

NAME _____

DAMPF BRAUN STEAM BROWN

REF 04116

E

STEAM BROWN

Wf

Name _____

Init. _____

Inhalt _____

Contents _____

LOT

REF 04112

F

REF 01116

STER. DAT. _____

STER. DATE _____

VERFALL DAT. _____

EXP. DATE _____

STER. NR. _____

STER. NO. _____

NAME _____

DAMPF BRAUN STEAM BROWN

STERIL

G

Steri Dat. _____

Steri Nr. _____

Name _____

Setname _____

Dampf Braun

REF 04003

H

INSERT LABELS

Insert labels make documentation easier for you, they can simply be pushed into the opening provided on the container. Available in different designs suitable for Aesculap, Martin, or other container systems. We also offer blank versions and with documentation fields in which sterilization data, batch numbers, and sieve names can be noted. Our insert labels are made of special cardboard compatible with sterilization. Single use.

DOCUMENTATION LABELS FOR WAGNER CONTAINERS

With indicator and documentation field.

A
B
C


REF	SELF-ADHESIVE	FORMAT	QTY
01133	no	40 x 59mm	1000 pcs
01134	yes	40/10 x 59mm	1000 pcs
01169	yes	40 x 59/74mm	1000 pcs

DOCUMENTATION LABELS FOR OTHER CONTAINERS

D
E
F
G
H

REF	DESCRIPTION	FORMAT	QTY
04067	with indicator, with documentation field, self-adhesive; for Geuder containers	23 x 60mm	1000 pcs
04116	with indicator, with documentation field, not self-adhesive; for CBM containers	39 x 55mm	1000 pcs
04112	with indicator, with documentation field, self-adhesive; for Microstop containers	39 x 63mm	1000 pcs
01116	with indicator, with documentation field, not self-adhesive; for sterilization containers	38 x 67mm	1000 pcs
04003	set label with indicator, self-adhesive, carrier material slotted from special cardboard 280gsm	40 x 80mm	1000 pcs

☆ Top selling items



DOCUMENTATION LABELS FOR AESCULAP CONTAINERS

For single use.

	REF	DESCRIPTION	FORMAT	QTY
☆ A	04102	blank, not self-adhesive	39 x 68mm	1000 pcs
B	01400	with indicator, with documentation field, not self-adhesive	35 x 80mm	1000 pcs
C	04184	without indicator, blank, not self-adhesive	35 x 80mm	1000 pcs
☆ D	01135	with indicator, with documentation field, not self-adhesive	39 x 68mm	1000 pcs
E	01136	with indicator, with documentation field, self-adhesive	39 x 67mm	1000 pcs
☆ F	04000	with indicator, with documentation field, self-adhesive	18 x 60mm	1000 pcs
G	04005	with indicator, with documentation field, not self-adhesive, 100 pieces each packed	18 x 50mm	1000 pcs
	04107	with indicator 'STERILE', without frame, self-adhesive	39 x 67mm	1000 pcs

DOCUMENTATION LABELS FOR MARTIN CONTAINERS

With steam indicator and documentation field. Single use.

	REF	SELF-ADHESIVE	FORMAT	QTY
H	01139	yes	46 x 50mm	1000 pcs
I	01122	yes	46 x 60mm	1000 pcs

DOCUMENTATION LABELS FOR MARTIN CONTAINERS

For single use.

☆

	REF	DESCRIPTION	FORMAT	QTY
J	01120	with indicator, with documentation field, self-adhesive	36 x 67mm	1000 pcs
K	04123	blank, not self-adhesive	36 x 67mm	1000 pcs
L	01121	with indicator, with documentation field, not self-adhesive	36 x 64 mm	1000 pcs

☆ Top selling items



A



B



C

TAG LABELS

Our tags help you to quickly label your containers and trolleys. They are heat-resistant up to 134 °C, not self-adhesive. Single use.



A

B

C

D

REF	DESCRIPTION	FORMAT	QTY
04162	with hole, white, made of polyester	90 x 50 x 20mm	1000 pcs
04192	with hole, red on one side, made of polyester	90 x 50 x 20mm	1000 pcs
04193	with hole, blue on one side, made of polyester	90 x 50 x 20mm	1000 pcs
04216	with hole, white, made of polyester	110 x 60 x 20mm	1000 pcs
04289	with hole, red on both sides, made of polyester	110 x 60 x 20mm	1000 pcs
04172	with cross-cut, white, made of polyester	75 x 40 x 20mm	1000 pcs



D

☆ Top selling items

Instrument Identification

The KeyDot® is a small, laser-etched, 2D Data Matrix barcode label that can be applied directly to an instrument's surface. Tested for sterilization efficacy and bio-compatibility for use on stainless steel surgical instruments, it has been conveniently incorporated into most of the instrument tracking software applications available today. With the KeyDot there is no specialized training or intricate equipment needed so application can be done on-site, eliminating the hassle of sending your instrument sets out for this tracking solution.

Our Identification Tape is designed to help you organize and identify instruments as they move between sterile processing and various departments. Using tape can help save you time, effort, and money by offering effective organization and identification of your surgical instrumentation. It is easy to apply and is compatible with steam and H₂O₂ sterilization.

We have the broadest selection of tape available and it is colorful, easy to use, and economical. Adhesive Remover is designed to dissolve tape residue on instruments, trays, basins, and other surfaces. It is the citrus based solution you need to remove residue without causing any damage.

Color Coding

Adhesive Remover

Roll Tape

Sheet Tape

KeyDot

Repair Tags


Printed Locking Tag

Cable Binders

Bundling Strings

Silicone Tubes






Color coding simplifies the assignment of individual products to different application or specialist areas or stations. We offer a wide range of identification options that can be selected according to preference or circumstances.

COLORED BEADS

Our colored beads made of silicone are ideal for marking instruments and are heat-resistant up to 134 °C. Available in various colors. Single use.


	REF	COLOR	MAX. Ø DIA	QTY
A	09652	turquoise	20mm	100 pcs
B	09653	green	20mm	100 pcs
C	09654	blue	20mm	100 pcs
D	09655	yellow	20mm	100 pcs
E	09656	black	20mm	100 pcs
F	09657	red	20mm	100 pcs
G	09658	white	20mm	100 pcs
H	09659	orange	20mm	100 pcs
I	09660	brown	20mm	100 pcs
J	09661	gray	20mm	100 pcs
K	09662	violet	20mm	100 pcs



COLORCURLS®

ColorCurls made of polypropylene are ideal for marking instruments and are heat resistant up to 134 °C. Available on request in various colors. For single use.

	REF	COLOR	RING HANDLE Ø DIA	QTY
A	09639	light blue	up to approx. 6mm	500 pcs
B	09640	dark blue	up to approx. 6mm	500 pcs
C	09641	light green	up to approx. 6mm	500 pcs
D	09642	dark green	up to approx. 6mm	500 pcs
E	09643	white	up to approx. 6mm	500 pcs
F	09645	dark red	up to approx. 6mm	500 pcs
G	09646	light yellow	up to approx. 6mm	500 pcs
H	09647	orange	up to approx. 6mm	500 pcs
I	09648	brown	up to approx. 6mm	500 pcs
J	09649	black	up to approx. 6mm	500 pcs
K	09651	turquoise	up to approx. 6mm	500 pcs



FIXING RINGS

Developed with an easy-to-use lock for marking sieve baskets and instruments or for attaching basket labels. These are made of plastic and are heat resistant up to 134 °C. For single use.

	REF	COLOR	INSIDE Ø DIA	QTY
A	09418	blue	10mm	100 pcs
B	09419	green	10mm	100 pcs
C	09420	white	10mm	100 pcs
D	09421	red	10mm	100 pcs
E	09422	yellow	10mm	100 pcs
F	09423	orange	10mm	100 pcs
G	09424	violet	10mm	100 pcs
H	09425	black	10mm	100 pcs
I	09426	transparent	10mm	100 pcs
J	09585	pink	10mm	100 pcs



ADHESIVE REMOVER

Our adhesive remover is designed to dissolve tape residue on instruments, trays, basins, or other surfaces. Biodegradable solution with citrus scent. Product availability and local country registration varies by market. Please contact your sales representative for information.

REF	DESCRIPTION	QTY
BC2006	adhesive remover 250ml bottle	6 bottles per case

KEY FEATURES:

- Removes adhesive residue from most surfaces
- 250ml bottle
- Citrus scent
- CE-marked
- pH 6 - 6.5
- Storage: place the item in a cool, dry area in its original container.
Keep container tightly closed when not in use
- Shelf life: 24 months from date of manufacture
- Transport hazard class 9

☆ Top selling items



ROLL TAPE

Available in different colors and widths. Compatible with steam sterilization up to 134 °C and H₂O₂. Single use.


	REF	COLOR	FORMAT	QTY
☆ A	T0102	blue	3.2 x 6350mm	1 roll
	T0202	blue	6.4 x 6350mm	1 roll
	T0302	blue	9.4 x 6350mm	1 roll
B	T0108	brown	3.2 x 6350mm	1 roll
	T0208	brown	6.4 x 6350mm	1 roll
	T0308	brown	9.4 x 6350mm	1 roll
☆ C	T0105	yellow	3.2 x 6350mm	1 roll
	T0205	yellow	6.4 x 6350mm	1 roll
	T0305	yellow	9.4 x 6350mm	1 roll
☆ D	T0103	green	3.2 x 6350mm	1 roll
	T0203	green	6.4 x 6350mm	1 roll
	T0303	green	9.4 x 6350mm	1 roll
E	T0107	orange	3.2 x 6350mm	1 roll
	T0207	orange	6.4 x 6350mm	1 roll
	T0307	orange	9.4 x 6350mm	1 roll
☆ F	T0104	red	3.2 x 6350mm	1 roll
	T0204	red	6.4 x 6350mm	1 roll
	T0304	red	9.4 x 6350mm	1 roll
G	T0109	black	3.2 x 6350mm	1 roll
	T0209	black	6.4 x 6350mm	1 roll
	T0309	black	9.4 x 6350mm	1 roll
H	T0106	violet	3.2 x 6350mm	1 roll
	T0206	violet	6.4 x 6350mm	1 roll
	T0306	violet	9.4 x 6350mm	1 roll
☆ I	T0101	white	3.2 x 6350mm	1 roll
	T0201	white	6.4 x 6350mm	1 roll
	T0301	white	9.4 x 6350mm	1 roll

☆ Top selling items


SOLID COLOR SHEET TAPE - 3.17MM WIDTH

Compatible with steam and H₂O₂ sterilization. Customizable.

Each 216mm x 279.4mm sheet includes:
462 pieces of 3.17mm x 25.4mm
88 pieces of 3.17mm x 38.1mm
22 pieces of 3.17mm x 50.8mm
22 pieces of 3.17mm x 127mm
(total of 1,899cm of tape)




	Number	Color
A	T11001	white
B	T11002	blue
C	T11003	green
D	T11004	red
E	T11005	yellow
F	T11006	purple
G	T11007	orange
H	T11008	brown
I	T11009	black
J	T11010	bright orange
K	T11011	pink
L	T11012	gray
M	T11013	lime green
N	T11014	light blue
O	T11015	lavender
P	T11016	beige
Q	T11017	plum




SOLID COLOR SHEET TAPE - 6.35MM WIDTH

Compatible with steam and H₂O₂ sterilization. Customizable.

Each 216mm x 279.4mm sheet includes:
231 pieces of 6.35mm x 25.4mm
44 pieces of 6.35mm x 38.1mm
11 pieces of 6.35mm x 50.8mm
11 pieces of 6.35mm x 127mm
(total of 950cm of tape)



	Number	Color
A	T12001	white
☆ B	T12002	blue
☆ C	T12003	green
☆ D	T12004	red
☆ E	T12005	yellow
F	T12006	purple
G	T12007	orange
H	T12008	brown
I	T12009	black
J	T12010	bright orange
K	T12011	pink
L	T12012	gray
M	T12013	lime green
N	T12014	light blue
O	T12015	lavender
P	T12016	beige
Q	T12017	plum
R	T12021	emerald green



☆ Top selling items


PATTERNED SHEET TAPE - 6.35MM WIDTH

Compatible with steam and H₂O₂ sterilization. Customizable. One 216mm x 279.4mm sheet pre-cut into 25.4mm, 38.1mm, 50.8mm and 127mm lengths. Cut as illustrated on page 194.

 T3200102 white/blue	 T3200103 white/green	 T3200104 white/red	 T3200105 white/yellow	 T3200106 white/purple
 T3200107 white/orange	 T3200108 white/brown	 T3200109 white/black	 T3200110 white/ bright orange	 T3200111 white/pink
 T3200113 white/ lime green	 T3200114 white/ light blue	 T3200115 white/lavender		
 T3200203 blue/green	 T3200204 blue/red	 T3200205 blue/yellow	 T3200206 blue/purple	 T3200207 blue/orange
 T3200211 blue/pink	 T3200503 yellow/green	 T3200504 yellow/red	 T3200506 yellow/purple	 T3200507 yellow/orange
 T3200511 yellow/pink	 T3200513 yellow/green	 T3200903 black/green	 T3200904 black/red	 T3200905 black/yellow
 T3200906 black/purple	 T3200907 black/orange	 T3200910 black/ bright orange	 T3200911 black/pink	 T3200912 black/gray

PATTERNED SHEET TAPE - 6.35MM WIDTH

Compatible with steam and H₂O₂ sterilization. Customizable. One 216mm x 279.4mm sheet pre-cut into 25.4mm, 38.1mm, 50.8mm and 127mm lengths. Cut as illustrated on page 194.

 T3200913 black/lime green	 T3200914 black/light blue	 T3200915 black/lavender	 T3200916 black/beige	
 T2200102 white/ blue squares	 T2200103 white/ green squares	 T2200104 white/ red squares	 T2200105 white/ yellow squares	 T2200106 white/ purple squares
 T2200107 white/ orange squares	 T2200109 white/ black squares	 T2200111 white/ pink squares	 T2200113 white/ lime green squares	
 T330010304 white/ green/red	 T330010902 white/ black/blue	 T330020503 blue/ yellow/green	 T330020504 blue/ yellow/red	 T330020703 blue/ orange/green
 T330030104 green/ white/red	 T330030504 green/ yellow/red	 T330040102 red/ white/blue	 T330040109 red/ white/black	 T330040203 red/ blue/green
 T330040309 red/ green/black	 T330040509 red/ yellow/black	 T330070103 orange/ white/green		

PATTERNED SHEET TAPE - 6.35MM WIDTH

Compatible with steam and H₂O₂ sterilization. Customizable. One 216mm x 279.4mm sheet pre-cut into 25.4mm, 38.1mm, 50.8mm and 127mm lengths. Cut as illustrated on page 194.

ER ER ER ER T50202 blue/ER	ER ER ER ER T50203 green/ER	ER ER ER ER T50204 red/ER	ER ER ER ER T50206 purple/ER	ER ER ER ER T50207 orange/ER
GI GI GI GI T5031609 beige/GI	OB OB OB OB T50411 pink/OB	OB OB OB OB T50414 light blue/OB	OR OR OR OR T52303 green/OR	OR OR OR OR T52307 orange/OR
DEN DEN DEN DEN T6030109 white/DEN	DEN DEN DEN DEN T60304 red/DEN	ENT ENT ENT ENT T60402 blue/ENT	ENT ENT ENT ENT T60403 green/ENT	ENT ENT ENT ENT T60404 red/ENT
ENT ENT ENT ENT T60407 orange/ENT	ENT ENT ENT ENT T60415 lavender/ENT	TX TX TX TX T50602 blue/TX	EYE EYE EYE EYE T60502 blue/EYE	EYE EYE EYE EYE T60503 green/EYE
EYE EYE EYE EYE T60504 red/EYE	EYE EYE EYE EYE T60506 purple/EYE	GEN GEN GEN GEN T60602 blue/GEN	GEN GEN GEN GEN T60603 green/GEN	GEN GEN GEN GEN T60604 red/GEN
GEN GEN GEN GEN T60606 purple/GEN	GEN GEN GEN GEN T60607 orange/GEN	GEN GEN GEN GEN T60609 black/GEN		

PATTERNED SHEET TAPE - 6.35MM WIDTH

Compatible with steam and H₂O₂ sterilization. Customizable. One 216mm x 279.4mm sheet pre-cut into 25.4mm, 38.1mm, 50.8mm and 127mm lengths. Cut as illustrated on page 194.

GYN GYN GYN GYN T60703 green/GYN	GYN GYN GYN GYN T60704 red/GYN	GYN GYN GYN GYN T60706 purple/GYN	GYN GYN GYN GYN T60707 orange/GYN	GYN GYN GYN GYN T60719 fuchsia/GYN
GYN GYN GYN GYN T60720 medical blue/GYN	ICU ICU ICU ICU T60803 green/ICU	ICU ICU ICU ICU T60804 red/ICU	L&D L&D L&D L&D T60903 green/L&D	L&D L&D L&D L&D T60904 red/L&D
L&D L&D L&D L&D T60911 pink/L&D	L&D L&D L&D L&D T60914 light blue/ L&D	NEU NEU NEU NEU T61002 blue/NEU	NEU NEU NEU NEU T61003 green/NEU	NEU NEU NEU NEU T61004 red/NEU
NEU NEU NEU NEU T61006 purple/NEU	NEU NEU NEU NEU T61007 orange/NEU	NEU NEU NEU NEU T61011 pink/NEU	PLA PLA PLA PLA T61203 green/PLA	PLA PLA PLA PLA T61204 red/PLA
SPD SPD SPD SPD T61403 green/SPD	SPD SPD SPD SPD T61407 orange/SPD	SPD SPD SPD SPD T61417 plum/SPD	PED PED PED PED T61114 light blue/PED	PED PED PED PED T61118 baby pink/PED
PED PED PED PED T6111309 lime green/ PED	PLA PLA PLA PLA T6120509 yellow/PLA	RAD RAD RAD RAD T6130709 orange/RAD	SPN SPN SPN SPN T61402 blue/SPN	

PATTERNED SHEET TAPE - 6.35MM WIDTH

Compatible with steam and H₂O₂ sterilization. Customizable. One 216mm x 279.4mm sheet pre-cut into 25.4mm, 38.1mm, 50.8mm and 127mm lengths. Cut as illustrated on page 194.

UROUROUROURO T61603 green/URO	UROUROUROURO T61604 red/URO	UROUROUROURO T61606 purple/URO	UROUROUROURO T61620 medical blue/URO	VASVASVASVAS T61704 red/VAS
VASVASVASVAS T61714 light blue/VAS	BARBARBARBAR T62308 brown/BAR	LOANLOANLOANLOAN T68603 green/LOAN	LOANLOANLOANLOAN T68607 orange/LOAN	CLNCCLNCCLNCCLNC T68720 medical blue/CLNC
ANESANESANESANES T68820 medical blue/ANES	NICUNICUNICUNICU T68903 green/NICU	NICUNICUNICUNICU T68907 orange/NICU		


PATTERNED SHEET TAPE - 6.35MM WIDTH

Compatible with steam and H₂O₂ sterilization. Customizable. One 216mm x 279.4mm sheet pre-cut into 25.4mm, 38.1mm, 50.8mm and 127mm lengths. Cut as illustrated on page 194.

♥♥♥♥ T7010104 white/hearts	♥♥♥♥ T70102 blue/hearts	♥♥♥♥ T70103 green/hearts	♥♥♥♥ T70104 red/hearts	♥♥♥♥ T70106 purple/hearts
♥♥♥♥ T70107 orange/hearts	♥♥♥♥ T70109 black/hearts	♥♥♥♥ T70111 pink/hearts	♥♥♥♥ T70121 emerald green/hearts	♥♥♥♥ T7010509 yellow/hearts
♥♥♥♥ T7011309 lime green/hearts	♥♥♥♥ T7011509 lavender/hearts	♥♥♥♥ T7020104 white/bones	♥♥♥♥ T70202 blue/bones	♥♥♥♥ T70203 green/bones
♥♥♥♥ T70204 red/bones	♥♥♥♥ T70206 purple/bones	♥♥♥♥ T70207 orange/bones	♥♥♥♥ T70208 brown/bones	♥♥♥♥ T70209 black/bones
♥♥♥♥ T70211 pink/bones	♥♥♥♥ T70219 fuchsia/bones	♥♥♥♥ T70221 emerald green/bones	♥♥♥♥ T7020509 yellow/bones	♥♥♥♥ T7021309 lime green/bones
♥♥♥♥ T7021509 lavender/bones	♥♥♥♥ T70302 blue/stars	♥♥♥♥ T70303 green/stars	♥♥♥♥ T70304 red/stars	♥♥♥♥ T70306 purple/stars
♥♥♥♥ T70307 orange/stars	♥♥♥♥ T70319 fuchsia/stars	♥♥♥♥ T70320 medical blue/stars	♥♥♥♥ T70321 emerald green/stars	




PATTERNED SHEET TAPE - 6.35MM WIDTH

Compatible with steam and H₂O₂ sterilization. Customizable. One 216mm x 279.4mm sheet pre-cut into 25.4mm, 38.1mm, 50.8mm and 127mm lengths. Cut as illustrated on page 194.

 T70406 purple/dots	 T70419 fuchsia/dots	 T70420 medical/ blue/dots	 T70421 emerald/ green/dots	 T70506 purple/triangles
 T70519 fuchsia/triangles	 T70520 medical/ blue/triangles	 T70521 emerald/g reen/triangles	 T70602 blue/ smiley faces	 T70603 green/ smiley faces
 T70604 red/ smiley faces	 T70606 purple/ smiley faces	 T7060509 yellow/ smiley faces	 T7060709 orange/ smiley faces	 T7061109 pink/ smiley faces
 T7061309 lime green/ smiley faces	 T7061409 light blue/ smiley faces	 T70702 blue/diamonds	 T70704 red/diamonds	 T70717 plum/diamonds
 T7080104 white/no	 T7090103 white/clovers	 T70903 green/clovers	 T7100118 Key Surgical logo	 T7110109 white/yin yang
 T7120104 white/lips	 T71302 blue/footprints	 T71303 green/footprints	 T71304 red/footprints	 T71306 purple/footprints
 T71307 orange/footprints	 T7131109 pink/footprints			

PATTERNED SHEET TAPE - 6.35MM WIDTH

Compatible with steam and H₂O₂ sterilization. Customizable. One 216mm x 279.4mm sheet pre-cut into 25.4mm, 38.1mm, 50.8mm and 127mm lengths. Cut as illustrated on page 194.

 T7140504 yellow/nuclear	 T7150109 white/ music notes	 T71603 green/crosses	 T71604 red/crosses	 T72002 blue/handprints
 T72004 red/handprints	 T7170103 white/pine trees	 T7180111 white/ribbon	 T71904 red/biohazard	
 T4010509 yellow/A	 T4020509 yellow/B	 T4030509 yellow/C	 T4040509 yellow/D	 T4050509 yellow/E
 T4070509 yellow/G	 T4120509 yellow/L	 T4140509 yellow/N	 T4160509 yellow/P	 T4161009 bright orange/P
 T4200509 yellow/T	 T4210509 yellow/U	 T8000509 yellow/0	 T8001009 bright orange/0	 T8010509 yellow/1
 T8020509 yellow/2	 T8030509 yellow/3	 T8040509 yellow/4	 T8050509 yellow/5	 T8060509 yellow/6
 T8070509 yellow/7	 T8080509 yellow/8			

PATTERNED SHEET TAPE - 6.35MM WIDTH

Compatible with steam and H₂O₂ sterilization. Customizable. One 216mm x 279.4mm sheet pre-cut into 25.4mm, 38.1mm, 50.8mm, and 127mm lengths. Each sheet contains numbers 1-44. Cut as illustrated on page 194.

1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22
1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22
1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22
1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22
23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44
23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44
23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44
23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44

T92002
blue/1-44


1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22
1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22
1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22
1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22
23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44
23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44
23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44
23	24	25	26	27	28	29	30	31	32	33	34	35	36	37	38	39	40	41	42	43	44

T92004
red/1-44

PATTERNED SHEET TAPE - 6.35MM WIDTH

Compatible with steam and H₂O₂ sterilization. Customizable. One 216mm x 279.4mm sheet pre-cut into 25.4mm and 38.1mm lengths. Cut as illustrated below.

Each 216mm x 279.4mm sheet includes:
308 pieces of 6.35mm x 25.4mm
44 pieces of 6.35mm x 38.1mm
(total of 950cm of tape)





T3100104
white/red



T3100105
white/yellow



T3100107
white/orange



T3100108
white/brown



T3100111
white/
light pink



T3100114
white/
light blue



T3100115
white/lavender



T3100502
yellow/blue



T3100503
yellow/
green



T3100901
white/black



T3100902
black/blue



T3100903
black/green



T3100904
black/red



T3100905
black/yellow



T3100910
black/orange



T3100911
black/light pink



T3100913
black/lime green



T3100917
black/plum



T330040902
red/
black/blue



T330040905
red/
black/yellow



T3200902
black/blue

PATTERNED SHEET TAPE - 6.35MM WIDTH

Compatible with steam and H₂O₂ sterilization. Customizable. One 216mm x 279.4mm sheet pre-cut into 25.4mm and 127mm lengths. Cut as illustrated on page194.

JAN
JAN
JAN
JAN

FEB
FEB
FEB
FEB

MAR
MAR
MAR
MAR

T90101
white/Jan
pink/Feb
light blue/Mar

APR
APR
APR
APR

MAY
MAY
MAY
MAY

JUN
JUN
JUN
JUN

T90201
green/Apr
orange/May
blue/Jun

JUL
JUL
JUL
JUL

AUG
AUG
AUG
AUG

SEP
SEP
SEP
SEP

T90301
red/Jul
brown/Aug
purple/Sep

OCT
OCT
OCT
OCT

NOV
NOV
NOV
NOV

DEC
DEC
DEC
DEC

T90401
lime green/Oct
yellow/Nov
lavender/Dec

A A A A A
B B B B B
C C C C C
D D D D D
E E E E E
F F F F F
G G G G G
H H H H H
I I I I I I
J J J J J J
K K K K K
L L L L L
M M M M M
N N N N N
O O O O O
P P P P P
Q Q Q Q Q
R R R R R
S S S S S
T T T T T
U U U U U
V V V V V
W W W W W
X X X X X
Y Y Y Y Y
Z Z Z Z Z

T91002
blue/alphabet

A A A A A
B B B B B
C C C C C
D D D D D
E E E E E
F F F F F
G G G G G
H H H H H
I I I I I I
J J J J J J
K K K K K
L L L L L
M M M M M
N N N N N
O O O O O
P P P P P
Q Q Q Q Q
R R R R R
S S S S S
T T T T T
U U U U U
V V V V V
W W W W W
X X X X X
Y Y Y Y Y
Z Z Z Z Z

T91004
red/alphabet

A A A A A
B B B B B
C C C C C
D D D D D
E E E E E
F F F F F
G G G G G
H H H H H
I I I I I I
J J J J J J
K K K K K
L L L L L
M M M M M
N N N N N
O O O O O
P P P P P
Q Q Q Q Q
R R R R R
S S S S S
T T T T T
U U U U U
V V V V V
W W W W W
X X X X X
Y Y Y Y Y
Z Z Z Z Z

T91019
fuchsia/alphabet

A A A A A
B B B B B
C C C C C
D D D D D
E E E E E
F F F F F
G G G G G
H H H H H
I I I I I I
J J J J J J
K K K K K
L L L L L
M M M M M
N N N N N
O O O O O
P P P P P
Q Q Q Q Q
R R R R R
S S S S S
T T T T T
U U U U U
V V V V V
W W W W W
X X X X X
Y Y Y Y Y
Z Z Z Z Z

T91003
green/alphabet

A A A A A
B B B B B
C C C C C
D D D D D
E E E E E
F F F F F
G G G G G
H H H H H
I I I I I I
J J J J J J
K K K K K
L L L L L
M M M M M
N N N N N
O O O O O
P P P P P
Q Q Q Q Q
R R R R R
S S S S S
T T T T T
U U U U U
V V V V V
W W W W W
X X X X X
Y Y Y Y Y
Z Z Z Z Z

T9100109
white/alphabet

PATTERNED SHEET TAPE - 6.35MM WIDTH

Compatible with steam and H₂O₂ sterilization. Customizable. One 216mm x 279.4mm sheet pre-cut into 25.4mm and 127mm lengths. Cut as illustrated on page 194.

AAAAAAA

T40119
red/A

♣♣♣♣♣
♥♥♥♥♥
♠♠♠♠♠
♦♦♦♦♦

T720040109
black/red/card suits

00000000
111111111
22222222
33333333
44444444
55555555
66666666
77777777
88888888

T9000109
white/0-8

00000000
111111111
22222222
33333333
44444444
55555555
66666666
77777777
88888888

T90003
green/0-8

00000000
111111111
22222222
33333333
44444444
55555555
66666666
77777777
88888888

T90006
purple/0-8

00000000
111111111
22222222
33333333
44444444
55555555
66666666
77777777
88888888

T90007
orange/0-8

00000000
111111111
22222222
33333333
44444444
55555555
66666666
77777777
88888888

T90019
fuchsia/0-8

00000000
111111111
22222222
33333333
44444444
55555555
66666666
77777777
88888888

T90020
medical blue/0-8

00000000
111111111
22222222
33333333
44444444
55555555
66666666
77777777
88888888

T90021
emerald green/0-8



HOW DOES KEYDOT WORK?

The KeyDot works like a license plate. It links a specific instrument to a data file within a software application, allowing that instrument to be tracked throughout the complete instrument processing cycle. Instruments marked with the KeyDot are scanned using a two-dimensional barcode reader. KeyDot works with many software applications available on the market today.

THE KEYDOT IS EFFICIENT

KeyDot provides you with the convenience to decide what to mark and when to mark it - you're in control. Since there are no equipment requirements or complicated procedures and training involved, instruments can be marked with a KeyDot on a continuous basis by anyone in your department, starting with critical instruments and expanding at an individual pace.

KEYDOT®

With the KeyDot there is no specialized training or intricate equipment required. Applying a KeyDot to a surgical instrument is quick and easy: simply clean, lift, place, and press. The instrument then needs to be sterilized to cure the adhesive and assure firm adhesion.


1 clean application area




2 lift KeyDot from sheet




3 place KeyDot on instrument



4 apply firm pressure and sterilize to secure adhesive





KEYDOTS SHOWN ACTUAL SIZE

KEYDOT®

	REF	Ø	QTY
☆ A	ID40501	3.2mm	1 sheet
	ID40502	3.2mm	2 sheets
	ID40505	3.2mm	5 sheets
	ID40510	3.2mm	10 sheets
☆ B	ID40001	4.8mm	1 sheet
	ID40002	4.8mm	2 sheets
	ID40005	4.8mm	5 sheets
	ID40010	4.8mm	10 sheets
☆ C	ID41001	6.4mm	1 sheet
	ID41002	6.4mm	2 sheets
	ID41005	6.4mm	5 sheets
	ID41010	6.4mm	10 sheets
☆ D	ID42001	9.5mm	1 sheet
	ID42002	9.5mm	2 sheets
	ID42005	9.5mm	5 sheets
	ID42010	9.5mm	10 sheets

CUSTOM NUMBERED KEYDOTS AVAILABLE ON REQUEST.



☆ Top selling items



A

B

C



D

E

F

LOCKING TAGS

Versatile, durable, and economical. Tamper-resistant. Compatible with steam and H₂O₂ sterilization. Single use. Customizable. 25.4mm length.

	REF	DETAILS	COLOR	QTY
A	LT150	CLEAN	white	100 pcs
	LT1501000	CLEAN	white	1000 pcs
B	LT210	QUICK TURNOVER	blue	100 pcs
	LT2101000	QUICK TURNOVER	blue	1000 pcs
C	LT310	LOANER	green	100 pcs
	LT3101000	LOANER	green	1000 pcs
★ D	LT450	REPAIR	red	100 pcs
	LT4501000	REPAIR	red	1000 pcs
E	LT550	SHARPEN	yellow	100 pcs
	LT5501000	SHARPEN	yellow	1000 pcs
F	LT720	BROKEN	orange	100 pcs
	LT7201000	BROKEN	orange	1000 pcs



★ Top selling items



REPAIR TAGS

These repair tags are available in plastic or Tyvek and are designed to mark defective instruments. We offer them blank or pre-printed (multilingual). Single use.

REF	DESCRIPTION	COLOR	IMPRINT	QTY
04132	made of Tyvek, washable up to 90 °C, 105 g/m ² , not suitable for steam sterilizer	white	blank	1000 pcs
04121	made of Tyvek, washable up to 90 °C, not suitable for steam sterilizer	white	blank	1000 pcs

Cable binders can be used to bundle different instruments. They can also be used to affix sieve basket labels or for quick visual identification of instruments, sets, and trays due to the unique properties. Our silicone products offer excellent fixation properties and at the same time protect the bundled objects. These facilitate the organization of the sieve during the packing process in the central sterile services department (CSSD) and during removal in the operating theatre. Available in different shapes, colors, lengths, and closure variants or different needs, such as in the operating theatre for marking instruments for repair or replacement. Tear-resistant. Single use.

CABLE BINDERS MADE OF SILICONE

Heat resistant up to 134 °C.

REF	COLOR	LENGTH	QTY
09246	white	60mm	100 pcs
09252	green	60mm	100 pcs
09253	blue	60mm	100 pcs
09254	yellow	60mm	100 pcs
09255	black	60mm	100 pcs
09256	red	60mm	100 pcs
09258	orange	60mm	100 pcs
09261	violet	60mm	100 pcs

A

B

C

D

E

F

G

H



A

B

C

D


E

F

G

H


☆ Top selling items



SILICONE CABLE BINDERS FOR BUNDLING INSTRUMENTS

Heat resistant up to 134 °C. Single use.

	REF	COLOR	LENGTH	QTY
A	09231	white	110mm	100 pcs
B	09235	red	110mm	100 pcs
C	09236	yellow	110mm	100 pcs
D	09237	green	110mm	100 pcs
E	09238	blue	110mm	100 pcs
F	09369	violet	110mm	100 pcs
G	09370	black	110mm	100 pcs
H	09373	white	110mm	100 pcs
I	09374	red	110mm	100 pcs
J	09375	yellow	110mm	100 pcs
K	09376	green	110mm	100 pcs
L	09377	blue	110mm	100 pcs
M	09378	orange	110mm	100 pcs
N	09379	black	110mm	100 pcs



CABLE BINDERS MADE OF SILICONE

Heat-resistant up to 134 °C. For single use.


	REF	COLOR	LENGTH	QTY
A	09282	white	190mm	100 pcs
B	09283	red	190mm	100 pcs
C	09284	green	190mm	100 pcs
D	09285	blue	190mm	100 pcs
E	09286	yellow	190mm	100 pcs
F	09381	orange	190mm	100 pcs
G	09382	black	190mm	100 pcs
☆ H	09239	red	210mm	100 pcs
☆ I	09240	yellow	210mm	100 pcs
☆ J	09241	green	210mm	100 pcs
K	09242	blue	210mm	100 pcs
L	09371	violet	210mm	100 pcs
M	09372	black	210mm	100 pcs
N	09586	white	210mm	100 pcs

PLASTIC CABLE BINDERS

Heat-resistant up to 134 °C. For single use.

	REF	COLOR	FORMAT	QTY
O	09482	white	200 x 3.5mm	100 pcs

☆ Top selling items



INSTRUMENT BUNDLING STRINGS

Made of silicone and heat-resistant up to 134 °C.

A



REF	FORM	COLOR	FORMAT	QTY
09143	flat	blue	300mm	100 pcs
09144	flat	white	300mm	100 pcs
09161	flat	yellow	500mm	100 pcs
09164	flat	gray	500mm	100 pcs
09039	flat	green	300mm	100 pcs
09139	flat	green	500mm	100 pcs
09146	flat	purple	300mm	100 pcs
09163	flat	red	500mm	100 pcs
09145	flat	black	300mm	100 pcs
09010	round	red-brown	Ø 3.0 x 700mm	100 pcs

TUBE BUNDLERS

B



REF	COLOR	Ø DIA	QTY
09158	white	25mm	100 pcs





A SILICONE TUBES

REF	DESCRIPTION	Ø INNER DIA	QTY
09043	transparent	1.0 x 0.3mm	25m/reel
09029	transparent	1.5 x 0.5mm	25m/reel
09037	transparent	1.5 x 0.8mm	25m/reel
09042	transparent	2.0 x 0.5mm	25m/reel
09028	transparent	2.0 x 0.8mm	25m/reel
09021	opaque red	2.0 x 1.0mm	25m/reel
09024	transparent	2.0 x 1.0mm	25m/reel
09149	transparent	2.5 x 1.0mm	25m/reel
09023	transparent	4.0 x 1.0mm	25m/reel
09115	transparent	4.0 x 1.0mm	25m/reel
09154	opaque blue	3.0 x 2.0mm	25m/reel
09022	transparent	5.0 x 1.0mm	25m/reel
09015	transparent green	5.0 x 1.0mm	25m/reel
09014	transparent	5.0 x 2.0mm	25m/reel
09041	transparent	6.0 x 3.0mm	25m/reel
09018	transparent	7.0 x 2.0mm	25m/reel



SILICONE TUBES

	REF	DESCRIPTION	Ø INNER DIA	CROSS-SECTION	QTY
A	09116	transparent	7.0 x 3.0mm	none	25m/reel
	09044	transparent	8.0 x 2.0mm	none	25m/reel
	09701	transparent	8.0 x 2.0mm	none	25m/reel
	09047	transparent	10 x 1.0mm	none	25m/reel
	09700	transparent	10 x 2.0mm	none	25m/reel
B	09710	transparent green, 6 ribs of 0.5mm each	2.5 x 0.5mm	6	25m/reel
	09706	transparent, 6 ribs of 0.5mm each	2.5 x 0.5mm	6	25m/reel
	09183	transparent, 6 ribs of 0.5mm each	2.5 x 1.2mm	6	25m/reel
	09233	transparent, 6 ribs of 1.0mm each	4.0 x 0.5mm	6	25m/reel
	09150	transparent, 6 ribs of 1.0mm each	4.0 x 1.2mm	6	25m/reel
	09151	transparent, 6 ribs of 1.0mm each	6.0 x 1.2mm	6	25m/reel
	09712	transparent green, 6 ribs of 1.0mm each	6.0 x 1.2mm	6	25m/reel
	09152	transparent, 6 ribs of 1.0mm each	8.0 x 1.2mm	6	25m/reel
	09142	transparent, 8 ribs of 1.0mm each	11 x 1.2mm	8	25m/reel
C	09153	transparent, 8 ribs of 1.0mm each	12 x 1.2mm	8	25m/reel
	09713	transparent green, 8 ribs of 1.0mm each	8.0 x 1.2mm	8	25m/reel
	09705	transparent, 16 ribs of 1.0mm each	20 x 1.0mm	16	25m/reel

QUALITY MONITORING AND ASSURANCE


Process monitoring and documentation have grown in importance in the field of medical device reprocessing over the past years which, in turn, has led to a significant increase in quality.

Routine controls in between annual maintenance and validation cycles provide consistent procedures and help to support compliance with the applicable guidelines and standards.

Monitoring Equipment for Sealing Devices and Accessories



A



INKCHECK® SEAL TIGHTNESS TEST

Suitable for routine testing of sealed seams in accordance with DIN EN ISO 11607-1, Annex B. Each set of InkCheck consists of 50 ampoules with 2 ml indicator solution and a magnifying glass card (REF 08529). InkCheck offers easy handling twist-off cap aiming to reduce spills and drips, suitable for all standard packaging paper/Tyvek®. The new generation of the seal seam tightness test InkCheck Plus provides simple and clean handling in one step. Single use.

A

B

REF	DESCRIPTION	QTY
08101	InkCheck	50 pcs
08116	InkCheck Plus	30 pcs
08124	InkCheck Plus Tyvek	30 pcs

SEAL PROOF

Seal Proof helps you validate sealers. Single use. Deviations in the sealing process, temperature and the contact pressure are shown on the indicator strip. All this in three simple steps:

- 1. Place Seal Proof in a cling-free clear pouch
- 2. Run through the sealing device once
- 3. Evaluate color change

C

REF	DESCRIPTION	FORMAT	QTY
08011	for transparent bags (paper foil)	173 x 91mm	250 pcs
08017	for Tyvek bags (Tyvek film)	172 x 90mm	100 pcs

B



C



DOCUMENTATION

Ever stricter standards mean that complete documentation is now part of everyday routine when reprocessing instruments. We offer a vast array of labels for identification of your sterile goods. Our range includes both paper labels for single use and longlasting foil-based labels. Take advantage of our comprehensive range and label your medical devices to meet your specific requirements.

Our foil and paper labels are available in almost all formats, and for paper labels you can choose between those with and without indicators for reliable identification in ongoing production for all common sterilization processes.

In our product range you will also find labels for low-temperature processes, such as H₂O₂ and thermo transfer foil in a selection of qualities.

Paper Labels

Foil Labels

Thermo Transfer Foil

Labelers and Labeling Accessories

Markers

Miscellaneous






A



B



C



D




E

Document the contents of your packaged medical devices using our high-quality paper labels. The labels are compatible with virtually all label printers, so you can print them according to your requirements and are thus always sure of finding the most important information at a glance. In our portfolio you will find labels in common formats. The duplex self-adhesive labels are particularly suitable for documentation in patient files. You will also find in our portfolio labels wound on the inside or outside and on various reel cores. Labels without tear-off perforations are ideal for cut-off printers.

SINGLE LABELS AND ENDLESS LABELS - WITHOUT INDICATOR

Duplex self-adhesive, blank without indicator, wound outside. For single use.


	REF	DESCRIPTION	TEAR-OFF PERFORATION	FORMAT	CORE- Ø DIA	QTY
A	04208	roll of 1000 labels	yes	50 x 33mm	40mm	1 roll
B	04084	suitable for non-woven materials, roll of 1000 labels	yes	70 x 33mm	40mm	1 roll
C	04152	roll of 1000 labels	yes	60 x 34mm	40mm	1 roll
D	04187	divided into 3 individual labels, roll of 2000 labels	yes	102 x 34mm 51 x 22mm 51 x 22mm	40mm	1 roll
E	04348	divided into 3 individual labels, roll of 500 labels	yes	102 x 34mm 51 x 22mm 51 x 22mm	40mm	1 roll
	04188	roll of 2000 labels	yes	75 x 35mm	40mm	1 roll
	04284	with slots, heat-resistant up to 134 °C, roll of 2000 labels	no	75 x 35mm	40mm	1 roll
	04161	roll of 2000 labels	no	63 x 40mm	40mm	1 roll




A

B


C



D



E



F



G

SINGLE LABELS - WITH STEAM INDICATOR (PINK TO BROWN)

Wound on the outside and duplex self-adhesive with tear-off perforation. Blank with steam indicator.

	REF	DESCRIPTION	FORMAT	CORE-Ø DIA	QTY
A	04010	with red colored corner - roll of 1000 labels	60 x 33mm	40mm	1 roll
B	04086	roll of 1000 labels	60 x 33mm	40mm	1 roll
C	04214	roll of 1000 labels	59 x 35mm	40mm	1 roll
D	04209	upper material with horizontal slits - roll of 1000 labels	59 x 35mm	40mm	1 roll
E	04290	roll of 1000 labels	58 x 38mm	40mm	1 roll
F	04141	roll of 1000 labels	57 x 40mm	40mm	1 roll
G	04167	with slit in the upper material - roll of 660 labels	72 x 48mm	40mm	1 roll
	04287	with slit in the upper material - roll of 660 labels	72 x 48mm	76mm	1 roll


☆ Top selling items



A



B



C



D



E



F

SINGLE LABELS - WITH STEAM INDICATOR SUITABLE FOR
NON-WOVEN MATERIAL

Blank with steam indicator (pink to brown). Wound outside and duplex self-adhesive.

	REF	DESCRIPTION	TEAR-OFF PERFORATION	FORMAT	CORE-Ø DIA	QTY
A	04160	roll of 1000 labels	yes	60 x 31mm	25mm	1 roll
B	04207	roll of 1000 labels	yes	50 x 33mm	40mm	1 roll
C	04089	roll of 1000 labels	no	70 x 33mm	40mm	1 roll
D	04295	roll of 1000 labels	no	70 x 33mm	76mm	1 roll
E	04301	with information on standards and color change roll of 1000 labels	yes	70 x 33mm	40mm	1 roll
	04336	with information on standards and color change, steam indicator left - roll of 1000 labels	yes	70 x 35mm	40mm	1 roll
F	04183	steam indicator right - roll of 1000 labels	yes	70 x 35mm	40mm	1 roll




☆ Top selling items



A



B



C



D



E



F

SINGLE LABELS - WITH STEAM INDICATOR SUITABLE FOR
NON-WOVEN MATERIAL

Blank with steam indicator (pink to brown). Wound outside and duplex self-adhesive.


- A
- B
- C
- D
- E
- F
- G

REF	DESCRIPTION	TEAR-OFF PERFORATION	FORMAT	CORE- Ø DIA	QTY
04250	roll of 1000 labels	yes	70 x 35mm	40mm	1 roll
04242	roll of 1000 labels	yes	70 x 35mm	76mm	1 roll
04273	with information on standards and color change, color bar red - roll of 1000 labels	yes	70 x 35mm	76mm	1 roll
04275	with information on standards and color change, green color bar - roll of 1000 labels	yes	70 x 35mm	76mm	1 roll
04276	with information on standards and color change, yellow color bar - roll of 1000 labels	yes	70 x 35mm	76mm	1 roll
04274	with information on standards and color change, blue color bar - roll of 1000 labels	yes	70 x 35mm	76mm	1 roll
04280	with information on standards and color change, vapor indicator left - roll of 1000 labels	yes	75 x 35m	76mm	1 roll
04241	roll of 1000 labels	yes	83 x 40mm	76mm	1 roll
04093	endless labels, wrapped outside, double self-adhesive, also sticks to non-woven materials	no	70 x 80mm	50mm	1 roll
04311	with information on standards and color change - roll of 1000 labels	yes	100 x 60mm	76mm	1 roll



G


☆ Top selling items



A



B



C



D



E



F

SINGLE LABELS - WITH STEAM INDICATOR (PINK TO BROWN)

Wound on the outside and single self-adhesive.

	REF	DESCRIPTION	FORMAT	CORE-Ø DIA	QTY
A	04190	suitable for non-woven material - roll of 1000 labels	60 x 31mm	40mm	1 roll
B	04094	suitable for non-woven material - roll of 2000 labels	60 x 31mm	40mm	1 roll
C	04265	roll of 1000 labels	60 x 33mm	40mm	1 roll
D	04198	for non-woven material - roll of 1000 labels	70 x 35mm	40mm	1 roll
E	04256	for non-woven material, indicator left - roll of 1000 labels	70 x 35mm	76mm	1 roll
F	04021	red, for non-woven material - roll of 2000 labels	90 x 35mm	40mm	1 roll



A




B




C



D



E



F

SINGLE LABELS - WITH STEAM INDICATOR (PINK TO BROWN)

Wound on the outside and single self-adhesive with tear-off perforation. With steam indicator.

REF	DESCRIPTION	FORMAT	CORE-Ø DIA	QTY
04299	for non-woven material, with information on standards - roll of 2000 labels	90 x 35mm	76mm	1 roll
04133L	for non-woven material and for longer warmth indicator left - roll of 750 labels	70 x 38mm	76mm	1 roll
04038	roll of 2000 labels	63 x 40mm	40mm	1 roll

ENDLESS LABELS - WITH STEAM INDICATOR (PINK TO BROWN)

Wound outside. Single use.

REF	DESCRIPTION	SELF-AD-HESIVE	ROLL LENGTH	WIDTH	CORE-Ø DIA	QTY
04085	suitable for non-woven material	duplex	80	70mm	50mm	1 roll
04145	suitable for non-woven material	duplex	80	70mm	76mm	1 roll
04142	with special adhesive (minimizes soiling of the cutting knife in the label printer)	duplex	80	70mm	50mm	1 roll
04268	with special adhesive	duplex	80	70mm	76mm	1 roll
04181	suitable for non-woven material, with free border on both sides	duplex	80	70mm	50mm	1 roll
04098	continuous label on reel with steam indicator	single	100	60mm	76mm	1 roll

☆ Top selling items



A



B



C



D



E



F

LABELS - WITH FO/EO AND STEAM INDICATOR

REF	DESCRIPTION	FORMAT	CORE-Ø DIA	QTY
04065	individual labels for Webeco, wound outside, perforated in the middle, with tear-off perforation, with back slit - roll of 1000 labels	86 x 36mm	40mm	1 roll
04045	continuous labels (EO = white, DI = pink), wound outside, double self-adhesive	70mm	50mm	1 roll

LABELS WITH H₂O₂ AND STEAM INDICATOR

With tear-off perforation, for single use, folded self-adhesive.

REF	DESCRIPTION	FORMAT	CORE-Ø DIA	QTY
04135	single labels, wound outside	70 x 35mm	40mm	1 roll
04279	single Tyvek labels, wound inside	70 x 35mm	76mm	1 roll
04097	continuous labels, wound outside	70mm	50mm	1 roll

OTHER SINGLE LABELS - WITH PUNCHED EDGE


With steam indicator (pink to brown), zigzag fold, with tear-off perforation. Single use.

REF	DESCRIPTION	FORMAT	QTY
04004	duplex self-adhesive, universal adhesive, printing: content, name, sterile date and patient.	99 x 37mm	2500 pcs

☆ Top selling items



A



B



C



D

Foil labels allow you to identify your containers and sieve trays. They offer durability and adhesion to smooth surfaces, even after the cleaning and disinfection process. We offer you different versions in numerous sizes compatible with the most common label printers.

FOIL LABELS - BLANK - WITHOUT INDICATOR

Heat-resistant up to 134 °C, suitable for washer-disinfectors, wound outside, single self-adhesive and tear-off perforation.

REF	DESCRIPTION	FORMAT	CORE-Ø DIA	QTY
04255	yellow, suitable for washer-disinfectors and dryers, strong adhesive	60 x 25mm	40mm	1 roll
04222	green, suitable for washer-disinfectors, strong adhesive	60 x 25mm	40mm	1 roll
04221	red, not washable, strong adhesive	60 x 25mm	40mm	1 roll
04227	white, suitable for washer-disinfectors, strong adhesive	50 x 23mm	40mm	1 roll
04176	white, suitable for washer-disinfectors, strong adhesive	60 x 25mm	40mm	1 roll
04186	white, suitable for washer-disinfectors, with one-sided barrel flap, can be removed without leaving any residue	55 x 35mm	50mm	1 roll



FOIL LABELS FOR IDENTIFYING AND SEALING CONTAINERS

It can be removed without any residue, single self-adhesive, wound outside. Single use.

REF	DESCRIPTION	FORMAT	CORE-Ø DIA	QTY
01001	with steam indicator	102 x 25mm	76mm	1 roll
01002	without steam indicator	102 x 25mm	76mm	1 roll
04096	without steam indicator	75 x 35mm	76mm	1 roll

☆ Top selling items

A



B

A THERMO TRANSFER FOILS

Thermal transfer reels with various coatings and features, suitable for your printer and labels. We recommend the premium resin which is particularly suitable for film carrier material (except 09858, which is only for paper). Wax/resin features are suitable for the most common requirements of the central sterile services department (CSSD) and especially for paper labels, and the wax feature is suitable for simple applications without thermal requirements. Coated on the outside or inside, so you are guaranteed to find the right foil reel for your printer in our range.

REF	COATING	QUALITY	FORMAT	SLEEVE LENGTH	SLEEVE-Ø DIA	STEAM STERILIZATION	QTY
09476	outside	resin	85 x 74mm	110mm	12.5mm	up to 150 °C	1 roll
09213	outside	wax	85 x 74mm	110mm	12.5mm	yes	1 roll
09452	outside	wax/resin	85 x 74mm	110mm	12.5mm	no	1 roll
09212	outside	wax	110 x 74mm	110mm	25mm	yes	1 roll
09853	outside	wax/resin	110 x 74mm	110mm	12.5mm	no	1 roll
09080	outside	wax	70 x 300mm	70mm	25mm	yes	1 roll
09270	outside	wax/resin	70 x 300mm	70mm	25mm	yes	1 roll
09069	inside	wax	76 x 300mm	76mm	25mm	yes	1 roll
09852	inside	wax/resin	76 x 300mm	76mm	25mm	no	1 roll
09486	outside	resin	80 x 300mm	80mm	25mm	up to 150 °C	1 roll
09854	outside	wax/resin	80 x 300mm	80mm	25mm	no	1 roll

ADHESIVE LABELS

Adhesive labels come in a DIN format. Single, self-adhesive.

B

REF	DESCRIPTION	INDICATOR	FORMAT	QTY
04016	for laser printers, single, self-adhesive - 25 x 44-label DIN A4 sheets	steam	48 x 25mm	1100pcs




A



B



C



D

SINGLE LABELS FOR LABELING DEVICES

REF	DESCRIPTION	COLOR	FORMAT	QTY
A 04072	METO - for labeling device 04019, with steam indicator	blue	29 x 28mm	12 rolls
04073	METO - for labeling device 04019, with steam indicator	yellow	29 x 28mm	12 rolls
04071	METO - for labeling device 04019, with steam indicator	green	29 x 28mm	12 rolls
04068	METO - for labeling device 04019, with steam indicator	orange	29 x 28mm	12 rolls
04070	METO - for labeling device 04019, with steam indicator	red	29 x 28mm	12 rolls
04069	METO - for labeling device 04019, with steam indicator	white	29 x 28mm	12 rolls
B 04054	for labeling device REF 04014, without indicator	blank	23 x 31mm	12 rolls
C 04041	for labeling device REF 04014, without indicator	blue	23 x 31mm	12 rolls
04040	for labeling device REF 04014, without indicator	green	23 x 31mm	12 rolls
D 04052	for labeling device REF 04014, with steam indicator	blue	23 x 31mm	12 rolls
04053	for labeling device REF 04014, with steam indicator	yellow	23 x 31mm	12 rolls
04048	for labeling device REF 04014, with steam indicator	orange	23 x 31mm	12 rolls
04050	for labeling device REF 04014, with steam indicator	red	23 x 31mm	12 rolls
04049	for labeling device REF 04014, without indicator	white	23 x 31mm	12 rolls



A



B



C



D



E



F

LABELERS AND ACCESSORIES

	REF	DESCRIPTION	QTY
A	04018	MONARCH - ink refill for "Monarch" labeler REF 04014	2 pcs
B	04119	METO ink refill for "Meto" labeler REF 04019	2 pcs
C	04014	MONARCH "Monarch" labeler - 3 lines	1 pc
D	04019	METO "Meto" labeler - 3 lines (12 characters per line (figures only); optionally available with text and special characters)	1 pc
E	04231	labeler for manual marking - 2 lines (10 characters per line -figures only)	1 pc
F	04230	labeler for manual marking - 3 lines (11 characters per line -figures only)	1 pc

A



B



C



D



E



F



G



H



I



J



K



L



MARKERS

	REF	DESCRIPTION	INK COLOR	LINE THICKNESS	QTY
A	05504	edding 8400 permanent marker, lightfast	black	0.5 - 1mm	10 pcs
	05513	permanent marker, sterilization-proof	red	0.75mm	10 pcs
B	05505	permanent marker, sterilization-proof	black	0.75mm	10 pcs
C	09089	edding 8010 X-ray marker	white	0.8mm	10 pcs
D	09087	edding 400 permanent marker, lightfast	blue	1mm	10 pcs
E	09170	edding 400 permanent marker, lightfast	red	1mm	10 pcs
F	09086	edding 400 permanent marker, lightfast	black	1mm	10 pcs
G	09091	edding 8020 skin marker, with round tip, not water resistant	blue	1mm	10 pcs
H	09090	edding 8020 skin marker, with round tip, not water resistant	red	1mm	10 pcs
I	05514	SteriLiner permanent marker, sterilization-proof	black	1mm	10 pcs
J	05515	SteriLiner XL permanent marker, with chisel tip, sterilization-proof	black	2 - 5mm	10 pcs
K	RT209P	Sharpie lab marker	black		12 pcs

BROTHER P-TOUCH

REF	MODEL	DESCRIPTION	QTY
L	09500	PTD600VP • German QWERTZ keyboard • color display, 3 lines x 20 characters • up to 24mm tape width • 280 displayable characters, phrase memory with 2800 characters • automatic saving of the last entry • 14 fonts, 10 print styles, 6 - 48 point size, 99 frame types and background patterns • vertical print, multiple printouts up to 99 times, multiline printouts up to 7 lines, mirror writing, tabulator function • 25 pre-programmed label templates (+ downloadable template), adjustable label length, 9 standard barcodes, multi-block printing up to 5 blocks • automatic cutter • recognizes the color of the inserted tape • USB Interface Layout design on the PC • logo and symbol printing in up to 180 x 360 dpi	1 pc

☆ Top selling items



A

B

C

D



E

F

G

H

LABELING TAPES FOR THE P-TOUCH LABELING MACHINE

Heat-resistant up to 134 °C, washable, only suitable for steam sterilization. Run length 8 meters.

	REF	TAPE WIDTH	LAMINATED	RIBBON COLOR	FONT COLOR	QTY
A	09511	18mm	yes	white	black	1 pc
B	09514	18mm	yes	yellow	black	1 pc
C	09544	18mm	no	white	black	1 pc
D	09543	24mm	no	white	black	1 pc
E	09502	24mm	yes	white	black	1 pc
F	09503	24mm	yes	red	black	1 pc
G	09505	24mm	yes	yellow	black	1 pc
H	09506	24mm	yes	green	black	1 pc

STENCIL TAPES FOR FOR CHEMICAL ETCHING

	REF	FONT COLOR	RIBBON COLOR	FORMAT	RUN LENGTH	QTY
I	09536	black	white	24mm	3m	1 pc

STATUS DOCUMENTATION

With tear-off perforation, for single use, double self-adhesive.

	REF	DESCRIPTION	FORMAT	QTY
J	09072	tray documentation block with steam indicator - block of 100 sheets	100 x 100mm	1 block

Achtung defektes Instrument/Gerät!
Containernummer: _____
Instrument: _____ Nr. _____
Art des Defektes
Instrument ☐ stumpf ☐ verbogen ☐ gebrochen
☐ schwergängig ☐ ausgeleiert ☐ abgenutzt
☐ Hartmetalleinlage beschädigt ☐ Oberfläche beschädigt
Maschine ☐ läuft nicht an ☐ läuft nicht rund ☐ undicht
☐ keine Kraft
Druckluftschlauch ☐ undicht ☐ Wandanschluß
☐ Maschinenanschluß
Sonstiges: _____
Datum: _____ Unterschrift: _____
Art.Nr. 09079 LOT

I

ZSVA
Name: _____ Datum: _____
fehlende Instrumente, z. ZL kein Ersatz:
dafür geliefert: _____
OP-SAAL
Name: _____ Datum: _____
Sieb-Nr.: _____ Kontakt-Sieb-Nr.: _____
Siebart: _____
fehlende Instrumente, Beanstandungen: _____
ZSVA
Name: _____ Datum: _____
fehlende Instrumente, nach Rücklauf: _____
REF: 09071

J

CSSD EQUIPMENT

On the following pages you will find various products for equipping your CSSD.

The importance of hospital logistics should not be underestimated when it comes to the reprocessing of medical devices and instruments. Short and long distances, safe transport and the avoidance of nosocomial (hospital-acquired) infections are all issues which need to be taken into account. Our UN certified transport boxes offer a safe system to protect personnel from used instruments during transport as well as sterile devices.


Dustproof bags are perfectly suited for transporting and storing sterile supplies. We provide dustproof bags in various sizes and material strengths, and with various sealing systems.

Dustproof Bags

Transport and Disposal

Cable Ties





DUSTPROOF BAGS WITH FLAP AND ADHESIVE SEAL

Solid design LDPE flat pouches. Transparent. Single use. Not suitable for sterilization.

☆ A

REF	DESCRIPTION	FORMAT	MATERIAL THICKNESS	QTY
03004	unprinted	100 x 150mm	50	500 pcs
03023	unprinted	150 x 270mm	50	100 pcs
03025	unprinted	150 x 380mm	50	100 pcs
03029	unprinted	200 x 300mm	50	100 pcs
03028	unprinted	200 x 650mm	80	100 pcs
03007	unprinted	300 x 400mm	60	100 pcs
03030	unprinted	300 x 500mm	50	100 pcs
☆ 03032	unprinted	600 x 750mm	50	100 pcs
03037	unprinted	630 x 820mm	60	100 pcs
03049	unprinted	700 x 1050mm	50	50 pcs
B 03033	unprinted	245 x 345mm	60	100 pcs
03036	unprinted	300 x 600mm	60	100 pcs
03017	unprinted	400 x 600mm	50	100 pcs
03031	unprinted	420 x 700mm	50	100 pcs
03022	unprinted	450 x 480mm	30	100 pcs
03035	unprinted	460 x 750mm	60	100 pcs
03034	unprinted	640 x 770mm	60	100 pcs
03040	with printing in German, English, Italian and French	500 x 500mm	60	100 pcs


DUSTPROOF BAGS WITH SLIDER, RESEALABLE - UNPRINTED

Solid design made of LDPE, flat pouch transparent. Single use. Not suitable for sterilization.

C

REF	FORMAT	MATERIAL THICKNESS	QTY
03041	320 x 230mm	75	100 pcs
03046	480 x 640mm	70	50 pcs
03045	320 x 440mm	70	100 pcs
03044	170 x 120mm	60	100 pcs
03043	200 x 150mm	60	100 pcs

☆ Top selling items



LDPE SECURITY BAG

For the hygienic storage after disinfection. The security closure is 40mm high and the tear-off label 50mm, with serial number with anti-tamper closure. For single use.

REF	FORMAT	INTERNAL DIMENSIONS	QTY
03077	520 x 500mm	500 x 500mm	100 pcs
03078	770 x 500mm	750 x 500mm	100 pcs

FLAT POUCHES / DUSTPROOF BAGS - UNPRINTED

For single use.

REF	DESCRIPTION	FORMAT	MATERIAL THICKNESS	QTY
09916	cleanroom quality, in secondary packaging	70 x 700mm	50	250 pcs





TRANSPORTATION BOXES

Stackable boxes that can be custom printed with text on the box and lid. They can withstand washer-disinfectors for decontamination purposes.

- BTB001
- suitable for transporting packing group I dangerous substances
- BTB002
- suitable for transporting packing group II dangerous substances
- BTB003
- suitable for transporting packing group I dangerous substances
- BTB004A
- suitable for transporting packing group II dangerous substances

	REF	WEIGHT	CAPACITY	COLOR	SIZE (EXTERNAL)	QTY
☆ A	BTB001	0.6kg	4l	blue	300 x 200 x 120mm	1 pc
B	BTB002	0.9kg	7l	blue	300 x 200 x 170mm	1 pc
☆ C	BTB003	1.4kg	16l	blue	400 x 300 x 220mm	1 pc
☆ D	BTB004A	3.3kg	49l	blue	600 x 400 x 320mm	1 pc
E	RTB002	0.9kg	7l	red	300 x 200 x 170mm	1 pc
☆ F	RTB003	1.4kg	16l	red	400 x 300 x 220mm	1 pc
☆ G	RTB004A	3.3kg	49l	red	600 x 400 x 320mm	1 pc

INNER BASKETS

Powder-coated aluminum inner baskets designed to fit transportation boxes.

REF	DESCRIPTION	QTY
BTB003BASK	inner basket to fit BTB003 and RTB003	1 pc
BTB004BASK	inner basket to fit BTB004A and RTB004A	1 pc

☆ Top selling items



TRANSPORTATION BOX OPENER

Stainless steel opener suitable for all BTB and RTB boxes.

REF	DESCRIPTION	QTY
BTBOPENER	box opener for BTB and RTB boxes, stainless steel	1 pc

TRANSPORTATION BOX TROLLEYS

Painted steel trolley with wheels for BTB and RTB boxes. Suitable for washer-disinfectors.

REF	DESCRIPTION	QTY
BTB004TROL	trolley to fit BTB004A and RTB004A	1 pc




TRANSPORTATION BOX LABELS

A selection of labels for 'BIOHAZARD' or 'USED MEDICAL DEVICE' to be inserted in plastic sleeves of transportation boxes. Our biohazard labels conform to the requirements for 6.2 hazard infectious substances.

REF	DESCRIPTION	SIZE	QTY
BTBLABEL003	biohazard, laminated label	125 x 40mm	20 pcs
BTBLABEL004	biohazard, laminated label	150 x 70mm	20 pcs
BTB232LRG	used medical device, laminated label	150 x 70mm	100 pcs





SECURITY CABLE TIES

Plastic cable ties in a range of colors for use during transportation and identification.

- A
- B
- ☆ C
- D
- E

REF	COLOR	SIZE	QTY
BTBREDTIE	red	3.6 x 145mm	100 pcs
BTBYELTIE	yellow	3.6 x 145mm	100 pcs
TIE001	blue	3.5 x 145mm	100 pcs
TIE006	clear	1.5 x 125mm	500 pcs
TIE010	red	4.8 x 300mm	100 pcs

CABLE TIES

Plastic cable ties in a range of colors for use during transportation. Colored cable ties help identification, along with numbered cable ties.

- F

REF	DESCRIPTION	COLOR	SIZE	QTY
TIE003	individually numbered	green	2 x 200mm	100 pcs
TIE005	individually numbered	red	2 x 200mm	100 pcs

TRANSPORTATION TRAYS AND STRAP

Transportation Trays reduce the risk of damage to wrapped packs due to poor handling and racking during transport and storage. Can be used with a strap or interlocked to form a lidded box.

- ☆ G
- H

REF	MATERIAL	SIZE	QTY
TTRAY002	polypropylene	500 x 270 x 70mm	1 pc
TTRAY003	cardboard	500 x 540 x 70mm	1 pc
STR004	blue strap	1400mm	1 pc

☆ Top selling items



Bestillinger sendes til:
Info@medonscandinavia.com